


## **Week 1: 8. June 2015 – 12. June 2015**

### **8. June 2015:**

#### 1) Amplification of *hxIA*, *hxIB* and *medh2* for TOPO cloning.

- Pipetting scheme and PCR program according to PCR with Phusion-HF DNA Polymerase Protocol
- Primer: *hxIA*\_P1/P2; *hxIB*\_P1/P2; *medh2*\_P1/P2
- Template: Genomic DNA from *Bacillus methanolicus* [100 ng/μl]



**Figure 1: Amplification of *hxIA*, *hxIB* and *medh2* for TOPO cloning.** 20 μl of PCR were checked on 1% (w/v) agarose gel. Expected sizes: *hxIA*-826 bp, *hxIB*-988 bp, *medh2*-1936 bp. As ladder (L) 1kb Ladder (NEB) was used.

#### 2) Liquid culture for growth of *Methylococcus capsulatus* (Bath)

- Inoculate a 30 ml culture of *Methylococcus capsulatus* Bath
- Use 22.5 ml NMS media and add 7.5 Methanol as a carbon source
- Incubate for 4 days at 37 °C

### **9. June 2015**

#### 1) Purification of generated *hxIA*, *hxIB* and *medh2* PCR products for TOPO cloning

- Using Wizard® SV Gel and PCR Clean-Up System (Promega) and following the provided manual.

#### 2) TOPO cloning of *hxIA*, *hxIB* and *medh2* into pCR4 Vector

- Using the Zero Blunt Topo cloning Kit (Invitrogen)
- See TOPO cloning protocol for pipetting scheme

#### 3) Transformation of *E. coli* TOP10 cells with TOPO cloning reactions


- Add 2 μl of TOPO reaction to chemically competent *E. coli* TOP10 cells.
- Following the protocol for transformation of chemical competent *E. coli* cells
- Plate on LB+Amp [100 μg/ml] and incubate overnight at 37 °C

### **10. June 2015**

- Transformation was successful and colonies were grown on the plate.

1) Colony-PCR to screen for clones containing the TOPO constructs

- Constructs: pCR4+*hxIA*, pCR4+*hxIB*, pCR4+*medh2*
- Pipetting scheme and PCR program according to PCR with Taq-DNA Polymerase Protocol
- Primer: *hxIA*\_P1/P2, *hxIB*\_P1/P2, *medh2*\_P1/P2
- Check 6 Clones per construct
- Positive Control: add 1 µl of PCR product of *hxIA*, *hxIB* and *medh2*
- Negative Control: add 1 µl MilliQ Water



**Figure 2: Colony-PCR to identify clones containing TOPO constructs.** Numbers 1-6 determine the checked clone. As positive control (PC) 1 µl PCR fragment was added. As negative control (NC) 1 µl MilliQ Water was added. 10 µl of PCR were analyzed on 1 % (w/v) agarose gel. Expected sizes: *hxIA*-826 bp, *hxIB*- 988 bp, *medh2*-1936 bp. As ladder (L) 1 kb Ladder was used.

2) Inoculate liquid culture for plasmid isolation of pCR4+*hxIA*, pCR4+*hxIB* and pCR4+*medh2*

- Inoculate clone 1 for plasmid isolation of pCR4+*hxIA* and clone 2 for plasmid isolation of pCR4+*medh2* with 5 ml LB+Amp [100 µg/ml] and incubate at 37 °C overnight shaking at 220 rpm.

- Although the Colony-PCR did not show positive pCR4+*hxIB* clones, we inoculated clone 3 with 5 ml LB+Amp [100 µg/ml], because also our positive control did not work, so we assumed a problem with the PCR reaction.


## **11. June 2015**

### 1) Plasmid Isolation of pCR4+*hxIA*, pCR4+*hxIB* and pCR4+*medh2*

- Use the QIAprep Spin Miniprep Kit (Qiagen) and follow the provided manual.

### 2) Amplification of *mmoC* and *mmoG* for TOPO cloning

- Pipetting scheme and PCR program according to PCR with Phusion-HF DNA Polymerase Protocol
- Primer: *mmoC*\_P1/P2, *mmoG*\_P1/P2
- Template: 1 µl of *Methylococcus capsulatus* liquid culture



**Figure 3: Amplification of *mmoC* and *mmoG* for TOPO cloning.** 20 µl of PCR were checked on 1% (w/v) agarose gel. Expected sizes: *mmoC*-1142 bp and *mmoG*-2249 bp. As ladder (L) 1kb Ladder (NEB) was used.

## **12. June 2015**

- Plasmids from iGEM2014 Team Braunschweig arrived. Each Plasmid consists of one gene encoding for the sMMO subunit (*mmoXYZBCD*) in the pSC1B3 vector backbone.

### Sequencing Reactions

- Send Plasmids pCR4+*hxIA*, pCR4+ *hxIB* and pCR4+*medh2* for sequencing.
  - Correct sequence of all cloned genes confirmed.
- Send Plasmids pSC1B3+ sMMO subunit gene for sequencing.
  - pSC1B3+*mmoY* contains a point mutation changing the codon GGT into GGC. Both codons encode the aminoacid Glycin.
  - Correct sequence of all other sMMO subunits confirmed.

