



**UNSW**  
AUSTRALIA

BABS UNSW iGEM Lab Protocol



Procedure	Name		Preparation of Agar Plates			
	Description		How to make Agar plates for bacterial growth			
Document	Name	Isabelle Capell-Hattam	Date	2/07/15	Version	1
Requirements	Time		~15 Minutes + Microwave and cooling			
	PPE		Gloves, Labcoat			
	Equipment		Microwave Pipettes and tips Sterile petri dishes Biosafety cabinet Parafilm Permanent Marker			
	Materials		Bacteriological growth media based agar (e.g LB agar for <i>E. Coli</i> , BG11 for <i>Synechocystis</i> ) Stock concentrations of antibiotics (that the bacteria have resistance to)			
Step 1	Melt the solid agar in a microwave till liquid and allow to cool until it can be held without burning					
Step 2	In the biosafety cabinet add 1 $\mu$ L of stock concentration antibiotic for every 1 mL of agar					
Step 3	Pour the agar into sterile petri dishes until $\sim\frac{1}{2}$ full. For a standard 10 cm dish this will be approximately 20 mL of agar.					
Step 4	Replace the cover on the dishes and allow to set at room temperature					
Step 5	When cool label the plates with the media and antibiotic they contain, seal individual plates in parafilm and place in a 4°C fridge until required for use					
Notes	Make sure the biosafety cabinet is sterile before pouring the plates					
Version History						