

## ***E.coli* Curdlan Purification**

### **Overview :**

Curdlan molecules may have as many as 12,000 glucose units and are insoluble in water, alcohols and most organic solvents, while it is soluble in dilute alkali (0.25 M NaOH).

We have two protocols for Curdlan purification from *E.coli* based on this properties to be soluble in NaOH solution.

### **Procedure n°1 :**

*Reference : Demonstration of Curdlan-type polysaccharide and some other in microorganisms with aniline blue, Itaru Nakanishi, Tokuya Harada and all, J. Gen. App!. Microbiol., 22,1-11(1976)*

1) The culture was mixed with an equal volume of 1 N NaOH and then centrifuged at 10,000 rpm for 10min to remove the cells.

2) The resulting supernatant was neutralized by adding 3 N HCl


3) The precipitate formed was collected by centrifugation at 10,000rpm for 10min.

-It was washed three times with water by centrifugation, and then dehydrated with acetone and dried in vacuo.

4) This precipitated polymer is the curdlan-type polysaccharide.

### **Procedure n°2 :**

*Reference : « Exopolymers from curdlan production: incorporation of glucose-related sugars by Agrobacterium sp. strain ATCC 31749. Jin W. Lee & all, Canadian Journal of Microbiology, 1997, Vol. 43, N° 2 : pages 149-156 »*



**Fig. 1.** Flow diagram depicting the purification procedure for curdlan. SN, supernatant. Adapted from the reference.

- 1) The culture was centrifuged at 12000 x g for 30 min.
- 2) The pellet was added to an equivalent volume of 0.5 N sodium hydroxide at 3°C
- 3) The mixture was stirred for 10 min
- 4) Left to stand for 3 h at the same temperature.
- 5) The resulting viscous solution was centrifuged at 12000 x g for 40 min
- 6) Curdlan in the clear supernatant was precipitated by neutralization with 10% acetic acid
- 7) Repeatedly washed with DWater, acetone and ether.