

**8th June, 2015**

- The following biobrick parts were revived from the **iGEM 2015 distribution plates** using the iGEM protocol Standard protocol.
- Transformation was carried with Ultra-competent cells for **DH5 $\alpha$**  strain of ***E.coli*** (with 200 $\mu$ L of LB).
- 200 $\mu$ L of the revived product was plated on LA+antibiotic.
- The plates were kept at 37°C overnight (12-16 hours).

Part Number	Location in Kit	Backbone
J23110 (Constitutive Promoter)	Plate 4, Well 19D	Ampicillin
K592025 (RBS+amilCP)	Plate 1, Well 19C	Chloramphenicol
J04450 (RFP)	Plate 4, Well 2H	Ampicillin

**Negative Control:** For both antibiotics (Amp and Cam)

**Positive Control:** RFP

**9th June, 2015**

Part Number	Number of colonies
J23110 (Promoter)	Many colonies
K592025 (RBS+amilCP)	Many single colonies
J04450 (RFP)	Lawn growth
Negative Control	No colonies

- Next, the colonies were inoculated in 5mL LB+ 0.1% antibiotic.
- They were kept in an incubator at 37°C overnight (12-16 hours) at around 180-200 rpm.

**10th June, 2015**

- Growth was observed in K592025, but no growth in J23110.
- Again a colony was picked from the plate with J23110 and inoculated in 5mL LB+ 0.1% antibiotic.
- They were kept in an incubator at 37°C overnight (12-16 hours) at around 180-200 rpm.

**11th June, 2015**

- Again no growth was observed in J23110 (even after inoculating twice).
- Transformation was done again for J23110
- Different promoter was also transformed.

Part Number	Location in Kit	Backbone
J23119 (Constitutive Promoter)	Plate 3, Well 17O	Chloramphenicol

### 13th June, 2015

- The two promoters were transformed with 200 $\mu$ L SOC media which was prepared using the iGEM protocol (earlier transformations with LB).
- 200 $\mu$ L of the revived product was plated on LA+antibiotic.
- The plates were kept at 37°C overnight (12-16 hours).

### 14th June, 2015

Part Number	Number of colonies
J23110 (Promoter 1)	30
K592025 (RBS+amilCP)	3
J23119 (Promoter 2)	47
Negative Control	No colonies

- The colonies obtained were inoculated in 5mL LB+ 0.1% antibiotic.
- They were kept in an incubator at 37°C overnight (12-16 hours) at around 180-200 rpm.

### 18th June, 2015

Plasmid isolation was carried out using the **Alkaline Lysis Protocol** from Sambrook and Maniatis.

Resuspended in TE Buffer.

### Results of Nanodrop:

Sample	Vial Number	A260/A280	Concentration (ng/ $\mu$ L)
K592025	1	1.99	6691.9
	2	2.11	4749.5
	3	2.12	4158.3
J23119	1	2.10	4865.3
	2	2.09	2757.4
	3	2.10	4902.3
J23110	1	2.09	5121.8
	2	2.15	8148.0
	3	2.13	3223.8

- The alkaline lysis buffer I wasn't added with RNase.
- The plasmid isolation was repeated with RNase added to Buffer I.

#### 27th June, 2015

- The colonies were again inoculated in 5mL LB+ 0.1% antibiotic.
- They were kept in an incubator at 37°C overnight (12-16 hours) at around 180-200 rpm.

#### 28th June, 2015

- RNase (10mg/mL stock) was added to Alkaline Lysis Buffer I.
- Resuspension in milliQ water.

#### Results of Nanodrop:

Sample	Vial Number	A260/A280	Concentration (ng/ $\mu$ L)
K592025	1	2.13	56.9
	2	1.39	67.4
	3	2.06	58.5
J23119	1	2.16	34.1

	2	2.15	32.3
	3	1.39	116.1
J23110	1	2.08	158.4
	2	2.02	3222.1
	3	1.95	580.4

### 29th June, 2015


- The colonies were again inoculated in 5mL LB+ 0.1% antibiotic.
- They were kept in an incubator at 37°C overnight (12-16 hours) at around 180-200 rpm.

### 30th June-10th July, 2015

Again, alkaline lysis was done.

### Results of Nanodrop:

Date	Sample	Vial Number	A260/A280	Concentration (ng/ $\mu$ L)
30/6	J23119	1	1.94	2376.8
30/6	J23119	2	1.95	2301.2
30/6	J23119	3	1.88	1662.9
8/7	J23119	1	1.99	1616.5
8/7	J23110	2	1.98	618.9
10/7	K592025	1	1.96	1641.1
10/7	K592025	2	1.97	1457.0
10/7	K592025	3	1.99	1404.6




- Legend:**
1. 2log ladder
  2. 3,17O (<2.1kb)-Promoter
  3. 3,17O (<2.1kb)-Promoter
  4. 3,17O (<2.1kb)-Promoter
  5. 3,17O (<2.1kb)-Promoter
  6. 4,19D (<2.1kb)-Promoter
  7. 1,19C (<2.8kb)-amilCP
  8. 1,19C (<2.8kb)-amilCP
  9. 1,19C (<2.8kb)-amilCP
  10. Blank



**26th August and 3rd September**


Results for Nanodrop by Qiagen spin mini kit

Sample	A260/A280	Concentration (ng/μL)
J23119	2.04	121.1
J23119	1.89	329.3
K592029	1.96	158
K592029	1.88	407.3



**Legend:**

- 1) ladder
- 2) 1,19C (chromoprotein)-cut (~2.9kb)
- 3) 1,19C (chromoprotein)-uncut (<2.9kb)
- 4) 3,170 (promoter)-cut(~2.1kb)
- 5) 3,170 (promoter)-uncut(<2.1kb)



**Legend:**

- 1) 2log DNA ladder
- 2) Other
- 3) Other
- 4) 1,19C (chromoprotein)-cut (~2.9kb)
- 5) 1,19C (chromoprotein)-uncut (<2.9kb)
- 6) 3,170 (promoter)-cut(~2.1kb)
- 7) 3,170 (promoter)-uncut(<2.1kb)

**26th August**

Single digestion reaction

(in $\mu$ L)	K592025	J23119
EcoRI	0.3	0.3
DNA	2	2
10XEcoRI Buffer	1	1
100X BSA	0.1	0.1
Distilled water	6.6	6.6
Total Volume	10	10

**27th August**

3A assembly double digestion

(in $\mu$ L)	Promoter(J23119)	amilCP(K592025)	Linear backbone
DNA	2	2	10
10X Buffer2	2.5	2.5	2.5
100X BSA	0.5	0.5	0.5
EcoRI	0.5	-	0.5
SpeI	0.5	-	-
PstI	-	0.5	0.5
XbaI	-	0.5	-
DpnI	-	-	0.5
Distilled water	14	14	5.5
Total Volume	20	20	20

Kept at 37°C for 4 hours, heat-kill at 80°C for 20min.

Colonies obtained (~10), were inoculated and miniprep was done with qiagen kit.

Table3: Double digestion for 2A assembly

2A assembly Double digestion:

(in $\mu\text{L}$ )	Promoter(J23119)	amilCP(K592025)
DNA	7.8	7
10X Buffer2	1	1
100X BSA	0.2	0.2
SpeI	0.5	0.5
PstI	0.5	-
XbaI	-	0.5
Distilled water	0	0.8
Total Volume	10	10

**9th September**