

iGEM 2015 – Microbiology – BMB – SDU

<p>Project type: Creation of a biobrick</p> <p>Project title: Linker-GFP</p> <p>Sub project: Assembly of a biobrick that contains the Linker-GFP gene</p>	<p>Creation date:</p> <p>Written by: ADK; EMT, CEM</p> <p>Performed by: ADK, EMT</p>
--	---

1. SOPs in use

iGEM2015_SOP010_v01_Phusion PCR
iGEM2015_SOP0014_v01_Gel purification
iGEM2015_SOP0017_v01_Fast digest
iGEM2015_SOP015_v02_ligation
iGEM2015_SOP009_v01_TSB transformation

2. Purpose

Assembly of a biobrick that contains the Linker-GFP gene.

3. Overview

Date (YY.MM.DD)	Person(s) (initials)	Experiments	SOPs
15.04.02	ADK, JSP	Phusion PCR x3	iGEM2015_SOP010_v01_P husion PCR
15.04.02		Gel purification	iGEM2015_SOP0014_v01_ Gel purification
15.04.03	ADK, JSP	Gel purification	iGEM2015_SOP0014_v01_ Gel purification
15.04.03	ADK, JSP, TBA, AC	Digest DNA-pieces	iGEM2015_SOP015_v02_ ligation
15.04.04	TBA, JSP, KBS	Transformation	iGEM2015_SOP009_v01_ TSB transformation

--	--	--	--	--

5. Experiment history

Date (YY.MM.DD)	Person(s) (initials)	Alterations to SOPs and remarks to experiments	SOPs
2015.04.02	ADK, JSP ADK, JSP JSP, TBA, AC	<p>First PCR attempt:</p> <ul style="list-style-type: none"> - Template: Linker_GFP - Primer: Linker-GFP_F + Con-mCherry_R - Optimised Temperature gradient: 55°C-64,9°C - Buffer: HF - Number of bp: about 1000 <p>Second PCR attempt:</p> <ul style="list-style-type: none"> - Template: Linker_GFP - Primer: Linker-GFP_F + Con-mCherry_R - Optimised Temperature gradient: 55°C-64,9°C - Buffer: 5 sample with HF and 4 sample with GC, - Number of bp: about 1000 <p>Gel purification:</p> <p>Used 15µl distilled water to elute the DNA</p> <p>Third PCR attempt:</p> <ul style="list-style-type: none"> - Template: Linker_GFP - Primer: Linker-GFP_F + Con-mCherry_R. Instead 2,5µl of each, we used 1,25µl of each - Optimised Temperature gradient: 55°C-64,9°C - Buffer: GC 	<p>-</p> <p>iGEM2015_SOP010_v01_P husion PCR</p> <p>-</p> <p>iGEM2015_SOP0014_v01_ Gel purification</p>

		- Number of bp: about 1000	
2015.04.03		<p>Gel purification: Used 15µl distilled water to elute the DNA</p> <p>Fast digest <u>PCR-product:</u> - Buffer: FastDigest - PCR product: [Linker-GFP]=24,8ng/µl, used 8µl - FastDigest enzyme: 1µl PstI + 1µl XbaI - Sterile water: 17µL</p> <p><u>Plasmid DNA:</u> Plasmid: pSB1C3, 114 ng/µl: 9 µl Reaction time: approx. 30 min</p> <p>Ligation of pieces of DNA: 3 samples: 1) <u>approx. 0 fmol Linker-GFP</u> - Plasmid (psB1C3 (R1)), 10 fmol/µl: 1µl - Linker-GFP: 0µl - H₂O: 16µl - DNA ligase buffer: 2µl - DNA ligase: 1µl 2) <u>approx. 20 fmol Linker-GFP</u> - Plasmid (psB1C3 (R1)), 10 fmol/µl: 1µl - Linker-GFP, 14 fmol/µl: 1,4µl - H₂O: 14,6µl - DNA ligase buffer: 2µl - DNA ligase: 1µl 3) <u>approx. 50 fmol Linker-GFP</u> - Plasmid (psB1C3 (R1)), 10 fmol/µl: 1µl - Linker-GFP, 14 fmol/µl: 3,6µl</p>	<p>iGEM2015_SOP0014_v01_ Gel purification</p> <p>iGEM2015_SOP0017_v01_ Fast digest</p> <p>iGEM2015_SOP0015_Ligati on of pieces of DNA</p>





		<ul style="list-style-type: none"> - H₂O: 12,4µl - DNA ligase buffer: 2µl - DNA ligase: 1µl <p>Leave the mixture overnight at 16°C</p>	
2015.04.04	TBA, KBC, JPS	<p>Transformation:</p> <p>Plasmid: linker-GFP with pSB1C3 backbone (0, 20, 50 fmol linker-GFP insert)</p> <p>E. coli wild type: MG1655</p> <p>Did not spin: approx. 100 µl was directly added on agar plate (Cml antibiotic)</p>	- iGEM2015_SOP009_v01_TSB transformation
2015.04.05	ADK, AC, KBC, JPS	<i>Same alterations as ligation made on 2015.04.03</i>	iGEM2015_SOP015_v02_ligation
2015.04.06	KBC	<p>Transformation:</p> <p>E. coli wild type: MG1655</p> <p>Plasmid: linker-GFP with pSB1C3 backbone (0, 20, 50 fmol linker-GFP insert): 15 µl added</p> <p>Approx. 100 µl was directly added on agar plate (Cml antibiotic)</p>	- iGEM2015_SOP009_v01_TSB transformation
2015.04.07	ADK	<p>Colony PCR with My Taq</p> <p>- 8 reaction in total, 5µl Taq, 2*1µl primer, 3µl steril water pr. reaction</p>	iGEM2013_SOP0021_v01 Colony PCR with my taq
15.04.11	TBA, AC	Fast Digest of linker-GFP from G8 cut with BamH1 and Pst1	iGEM2015_SOP0017_v01_Fast digest
15.04.11	TBA, AC	Gel purification , reveals that Fast digest of linker-GFP was unsuccessful.	iGEM2014_SOP0014_v01_Gelpurification
15.04.11	TBA, AC	<p>Fast digest (2nd attempt) cut with BamH1 and Pst1</p> <p>Gel purification (2nd attempt)</p>	iGEM2015_SOP0017_v01_Fast digest

		Successful Fast digest. Phusion PCR from G7, using UF-buffer and half amount of primers. temp gradient from 55 degC to 65 degC and left overnight	iGEM2015_SOP010_v01_P husion PCR
15.04.12	AC	Gel purification	iGEM2014_SOP0014_v01_ Gelpurification
15.04.19 15.04.19	TBA, AC, JSP	Fast digest x3 11,5 µL of G9 for total volume of 30 µL. Using BamH1 and Pst1, incubating for 30 min @37 degC runing a gel for 35 min Gel purification x3 eluted in 50 µL water, stored as Y12 = linker-GFP; 7,3 ng/µL Y13 = linker-GFP; 8,7 ng/µL Y14 = linker-GFP; 11,2 ng/µL	iGEM2015_SOP0017_v01_ Fast digest


15-05-06	AC	Ligation GFP with T18 and T25. See two hybrid protocol for further information	iGEM2015_SOP015_v02_ ligation

6. Results

Date (YY.MM.DD)	Picture	Comments


<p>2015.04.02</p>		<p>First PCR attempt. All we see are primer dimers</p>
		<p>Second PCR attempt. We see a lot of primer dimer and very weak 3 bonds about 1000bp long in the 3 last bonds with GC buffer. We purify the sample from the gel, and gets a concentration at 3,1ng/μl. We runs a PCR at the sample overnight, in aim to, increase the concentration, but this time with half as much primer, and only with GC buffer.</p>
<p>2015.04.03</p>		<p>Third PCR attempt. Successful attempt, we see 5 bonds about 1000bp long. We purified the sample, and got the concentration 24,8ng/μl.</p>
		<p>Digest DNA-pieces: (As wanted) is the GFP linker about 1000 bp long. We purified the sample and got the concentration 9,2ng/μl. The digest of the plasmid is successful. The bond about 1000bp is “waste”, The bond about 3000 is our successful digested plasmid backbone.</p>

2015.04.05




Transformation:
No colony was seen

2015.04.06




Transformation (from yesterday):
Red and white colonies at all agar plates. The white colonies at the 0-sample are religations.




Colony PCR.
The first bond didn't work (colony 1). The second bond (colony 2) was under 500bp, which shows it is a religation. The bonds 3-8 shows about 1200bp, which fit the length of the plasmid backbone and GFP-linker. TJECK the length of GFP linker!

15.04.11



Gel-results from 2nd Fast Digest, Succesful

15.04.19



Gel-results from Fast digest of overnight PCR

7. Appendices