

## Protocol for PCR Construction of HPV detector plasmid

1. HPV Detector Plasmid (uncut) linearised with KpnI for 4 hours at 37 degrees & heat inactivated at 80 degrees for 5 minutes

Sample	Volume (ul)
DNA (68ng/ul)	1.47 (100ng)
KpnI	0.5
Buffer 4	2
H2O	16.03
<b>Total</b>	<b>20</b>

2. Digestion was diluted 1 in 1000 ul of TE buffer (5pg DNA per ul)
3. PCR reaction was set up as follows:

### PCR Reactions :

	1	2	3 (Control)
DNA (5pg/ul)	4	4	4
Solis Biodyne 5x Hot Start	4	4	4
F Primer	1	1	0
R Primer	1	1	0
H2O	10	10	12
<b>Total</b>	<b>20</b>	<b>20</b>	<b>20</b>

### PCR Conditions:

Step	Temperature (degrees celsius)	Time	Number of
------	-------------------------------	------	-----------

			<b>Cycles</b>
<b>Initial denaturation</b>	95	15 mins	1
<b>Denaturation</b>	95	30 sec	35
<b>Annealing</b>	56	40 sec	35
<b>Elongation</b>	72	3 min	35
<b>Final Elongation</b>	72	7 min	1

1. Run 3ul of PCR products on gel to check if present
2. DpnI digest carried out by adding 0.5ul of DpnI to reaction tube with PCR product - incubate for 1 hour @ 37 degrees & 80 degrees for 10 mins
3. PCR clean up carried out after digest
4. Run 5 ul of clean up product on gel
5. Quadruple digest carried out on sample by adding 1 ul of NtBsQI + 5ul Buffer 4 and incubating for 2 hours at 50 degrees. The final 3 enzymes (KpnI, HindIII & NtBsl) were then added and incubation continued for 6 hours at 37 degrees
6. Digested reaction was reacted with decoy oligos by adding 1ul of 100mM decoys (R & L) and allowing to sit for 10 mins at room temperature
7. PCR Cleanup was carried out again.
8. Run 5ul of product on gel to check DNA is present

## **Testing PCR Constructed Detector**

### **Mixes Yielded**

	<b>Sample</b>	<b>Reaction</b>
<b>1</b>	PCR product with DpnI digest	A
<b>2</b>	PCR product without DpnI digest	
<b>3</b>	PCR Control (no detector)	
<b>4</b>	PCR product with DpnI digest, PCR cleanup, 4x digest	B & C
<b>5</b>	PCR product with DpnI digest, PCR cleanup, 4x digest, Decoy reaction, PCR cleanup	
<b>6</b>	Digested detector (made 10/6/15)	

7	Digested detector - working (made Nov. 2014)	
---	--	--


**16/6/15 - Reactions carried out prior to transformation:**

<b>T</b>	Target - Detector Reaction with Target
<b>C</b>	Control - Detector Reaction without Target
<b>Detector Reaction</b>	<p><b>Addition of GFP to Detector mixes 4 &amp; 5</b> 0.3 ul GFP added per 2ul of detector</p> <p><b>Detector Reaction</b> See other protocol 1ul of Detector + GFP used for each 1 ul of ds target used</p>

**Results**

- Results of transformation after 5ul of each reaction transformed on 16/6/15

Reaction	Non - GFP (Target)	GFP Expressing CFU
<b>6T</b>	16	38
<b>6C</b>	0	4
<b>4T</b>	0	1
<b>4C</b>	0	6
<b>5T</b>	13	68
<b>5C</b>	0	33



Transformation of cells with Digested Detector (right) and PCR Amplified detector (left). Plates with target on top, no target on bottom. GFP also added to reaction. (16/6/15)


**17/6/15 - Repeat of Transformation using Highly Competent Cells and Detectors 5 & 6**

- Detector reaction carried out as before
- Transformation followed NEB protocol for highly competent cells

**Results of transformation using highly competent cells and detectors 5 & 6 (PCR Mix 5 & Digested detector)**


	<b>NGFP</b>	<b>GFP</b>	<b>Ratio</b>
<b>Digested detector + target</b>	2,200	1,950	1.128
<b>PCR 5 + target</b>	761	838	0.908

<b>Digested Detector + no target</b>	11	1,500	0.007
<b>PCR 5 + no target</b>	0	1,500	0




Transformation of highly competent cells with digested detector (right) and PCR amplified detector (left) after reaction with target. Blue points note NGFP containing CFUs (not detector plasmid). (18/6/15)


### Testing PCR



Transformation of competent cells with 30bp HPV Target Detector (left) and 55bp HPV Target Detector (right). Blue points note target - positive CFUs. (22/6/15)



Transformation of competent cells with PCR Amplified SRY Target Detector. Control plate (left) shows 78 GFP expressing CFUs. Target plate (right) shows 89 NGFP expressing CFUs (detector + target) and 7 GFP CFUs. (22/6/16)



Transformation of competent cells with PCR Amplified 30bp HPV Target Detector. Target plate (left) shows 134 GFP expressing colonies. Control plate (right) shows 0 colonies. (22/6/15)

**Results of Transformation using PCR detectors (22/6/15)**

	GFP	NGFP
HPV PCR Amplified	8	0
HPV PCR Amplified + Target	188	12
HPV-GFP PCR Amp	0	0

<b>HPV-GFP PCR Amp + Target</b>	356	0
<b>HPV-GFP 30bp Target</b>	0	0
<b>HPV-GFP 30bp Target + Target</b>	~550	0
<b>SRY PCR Amp</b>	80	4
<b>SRY PCR Amp + Target</b>	14	178

24/6/15 Transformation of SRY Detector with Genomic DNA


<b>Tube</b>	<b>1 - Male (digest)</b>	<b>2 - Male (no digest)</b>	<b>3 - Female (digest)</b>	<b>4 - No DNA (control)</b>
<b>SRY Detector</b>	1	1	1	1
<b>GFP (0.01ul plasmid /ul)</b>	1	1	1	1
<b>Tris</b>	1	1	1	1
<b>Male - Digested DNA</b>	1	0	0	0
<b>Male - Undigested DNA</b>	0	1	0	0
<b>Female - Digested DNA</b>	0	0	1	0
<b>H2O</b>	6	6	6	7
<b>Total</b>	10	10	10	10



### Results of Transformation (25/6/15) & (26/6/15)

Tube	1 - Male (digest)			2 - Male (no digest)			3 - Female (digest)			4 - No DNA (control)		
NGFP	2	9	15	2	4	17	7	12	7	2	2	
GFP	17	1216	1064	15	856	1840	25	1128	2104	0	1572	
Adjusted Ratio	0.0059	0.007	0.014	0.0067	0.005	0.009	0.014	0.011	0.003	0	0.001	
Average	0.00897			0.0069			0.0093			0.0003		
Average 2	0.0105			0.007			0.007			0.0005		

Note : Test 1 used less Control GFP, Average 2 adjusts for this.



Transformation carried out on 25/6/15 using SRY Detector and genomic samples from male and female mice. Male digested DNA (Left), Male undigested DNA (Middle), Female digested DNA (Right). All samples tested in duplicate.