

iGEM 2015 – Microbiology – BMB – SDU

Project type: Screening		Creation date:
Project title: Peptide aptamer screening		Written by: TBA, EMT, CEM
		Performed by: TBA, JSP, AC and EMT

1. SOPs in use

iGEM2014_SOP0017_v01_ Fast digest

iGEM2014_SOP015_v02_ligation

iGEM2013_SOP0009_v01_TSB transformation

iGEM2013_SOP0021_v01_ Colony PCR with MyTaq

iGEM 2014 SOP0019 - Miniprep

2. Purpose

To screen for peptide aptamers using the two-hybrid system.

3. Overview

Date (DD.MM.YY)	Person(s) (initials)	Experiments	SOPs
11.08.15	TBA, JSP	Fast Digest of pSB1C3 and pSB1K3 with Eco + Pst and Xba + Pst. Fast Digest of “Scaffold-3xFLAG” with Eco + Pst Fast Digest of “Intein” with Xba + Pst	iGEM2014_SOP0017_v01_ Fast digest

		Ligation of Scaffold and Intein with the two backbone. Ligation left overnight	iGEM2014_SOP015_v02_ligation
12.08.15	TBA	Transformations of Ligations into MG1655	iGEM2013_SOP0009_v01_TSB transformation
13.08.15	TBA, JSP	Redone transformation from ligation-leftovers	iGEM2013_SOP0009_v01_TSB transformation
14.08.15	JSP	New ligation of Scaffold and Intein in pSB1C3 + pSB1K3	iGEM2014_SOP015_v02_ligation
15.08.15	ADK	Transformation of ligations	iGEM2013_SOP0009_v01_TSB transformation
16.08.15	TBA	Colony-PCR	iGEM2013_SOP0021_v01_Colony PCR with MyTaq
17.08.2015	ADK, TBA	Miniprep on overnight cultures Freeze stock preparation	Bio-Rad Quantum Prep: plasmid miniprep Kit
17.08.15	TBA, JSP	Fast digest of received targets proteins and appropriate backbones.	iGEM2014_SOP0017_v01_Fast digest
17.08.15	ADK, TBA	Ligations of target proteins, linker-intein-scaffold (LIS) and into different backbones	iGEM2014_SOP015_v02_ligation
18.08.15	ADK	Transformation of ligations from yesterday	iGEM2013_SOP0009_v01_TSB transformation
24.08.15	JSP	Ligation of PstA-RFP(Y79) and T25-IL2 (Y80) and T25-Hfq (Y83).	iGEM2014_SOP015_v02_ligation

25.08.15	JSP	Transformation of ligations from 24.08.2015.	iGEM2013_SOP0009_v01_TSB transformation
26.08.15	JSP	Colony-PCR on transformations from 25.08.2015	iGEM2013_SOP0021_v01_Colony PCR with MyTaq
27.08.15	JSP	Miniprep on PctA-RFP-T25-Hfq	Bio-Rad Quantum Prep: plasmid miniprep Kit
29.08.15	AC TBA	Fast Digest of pSB1C3-T18-LIS. Cut with EcoI + PstI. Stored as Y105 Ligation of digestion into pSB1K3-backbone	iGEM2014_SOP0017_v01_Fast digest iGEM2014_SOP015_v02_ligation
30-08.2015	EMT	Transformation	iGEM2013_SOP0009_v01_TSB transformation
31.08.2015	EMT	Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI + PstI. Ligation of digestion into pSB1K3-backbone Transformation	iGEM2014_SOP0017_v01_Fast digest iGEM2014_SOP015_v02_ligation iGEM2013_SOP0009_v01_TSB transformation
01.09.2015	EMT	Ligation of digestion from 31.08.2015 into pSB1K3-backbone Transformation	iGEM2014_SOP015_v02_ligation iGEM2013_SOP0009_v01_TSB transformation

2.09.2015	EMT	Kolony PCR on pSB1K3-T18-LIS pSB1K3-LSx3FLAG	iGEM2013_SOP0021_v0 1_Colony PCR with MyTaq
3.9.2015	EMT	Miniprep of overnight Cultures	iGEM 2014 SOP0019 - Miniprep
5.08.2015	EMT	Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with XhoI. Fast Digest of Library with XhoI Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library Transformation	iGEM2014_SOP0017_v01_ Fast digest iGEM2014_SOP015_v02 _ligation iGEM2013_SOP0009_v0 1_TSB transformation
06.08.2015	EMT	Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library Transformation	iGEM2014_SOP015_v02 _ligation iGEM2013_SOP0009_v0 1_TSB transformation
07.08.2015	EMT	Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with XhoI. Fast Digest of Library with XhoI Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library Transformation	iGEM2014_SOP0017_v01_ Fast digest iGEM2014_SOP015_v02 _ligation iGEM2013_SOP0009_v0 1_TSB transformation

08.09.2015	EMT	<p>Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI and PstI.</p> <p>Fast Digest of pSB1A2 with EcoI and PstI</p> <p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with pSB1A2</p> <p>Transformation</p>	<p>iGEM2014_SOP0017_v01_Fast digest</p> <p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>
09.09.2015	EMT	<p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with pSB1A2</p> <p>Transformation</p>	<p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>
10.09.2015	EMT	<p>Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI and PstI.</p> <p>Fast Digest of pSB1A2 with EcoI and PstI</p> <p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with pSB1A2</p> <p>Transformation</p> <p>Koloni PCR on transformation from 09.09.2015</p>	<p>iGEM2014_SOP0017_v01_Fast digest</p> <p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p> <p>iGEM2013_SOP0021_v01_Colony PCR with MyTaq</p>

11.09.2015	EMT	Koloni PCR on transformation from 10.09.2015	iGEM2013_SOP0021_v01_Colony PCR with MyTaq
12.09.2015	EMT, AC	<p>Restriction Blot, done on pSB1A2-T18-LSx3FLAG, cuts made;</p> <p>-EcorI -EcorI+PstI -XhoI and uncut pSB1A2</p> <p>Miniprep on T18-LSx3FLAG</p> <p>Fast Digest on pSB1A2-T18-LSx3FLAG with XhoI</p> <p>Ligation of digestion with Library</p> <p>Transformation of targets (T25-linker-CRSTA) into BTH101</p>	<p>iGEM 2014 SOP0019 - Miniprep</p> <p>iGEM2014_SOP0017_v01_Fast digest</p> <p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>

4. Materials required.

Materials in use

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations

5. Experiment history

Date (DD.MM.YY)	Person(s) (initials)	Alterations to SOPs and remarks to experiments	SOPs
11.08.2015	TBA, JSP	<p>Fast digest of Scaffold (using 15 μL of IDT-stock) with Eco + Pst. Stored as Y68</p> <p>Fast Digest of Intein (using 15 μL of IDT-stock) with Xba + Pst. Stored as Y69</p> <p>Fast digest of pSB1K3 (R71 using 10 μL) with Eco + Pst and Xba + Pst. Stored as Y70 and Y71, respectively.</p> <p>Fast digest of pSB1C3 (R72 using 5 μL) with Eco + Pst and Xba + Pst. Stored as Y72 and Y73 respectively.</p> <p>Backbone digestions run on gel, Scaffold and Intein purified using coloumbs.</p>	<p>iGEM2014_SOP0017_v01_ Fast digest</p> <p>iGEM2014_SOP015_v02 _ligation</p>

		<p>Otherwise done according to SOP.</p> <p><u>Ligation:</u> <u>Backbones to obtain 10 fmol</u> Y70: 1,6 µL was used Y71: 1,8 µL was used Y72: 1,33 µL was used Y73: 1,15 µL was used</p> <p><u>Scaffold (Y68)</u>, for 20 fmol 1,4 µL was used. for the 50 fmol ligation 3,5 µL was used.</p> <p><u>Intein (Y69)</u>, for 20 fmol 1,5 µL was used. for 50 fmol 4,5 µL was used</p> <p>Ligations left for ~20 hours</p>	
12.08.2015	TBA	Done according to SOP.	iGEM2013_SOP0009_v01_TSB transformation
13.08.2015	TBA, JSP	Previous transformation failed, possibly due to mess up during plating. Transformation redone	iGEM2013_SOP0009_v01_TSB transformation
14.08.2015	JSP	<p>Ligation of scaffold (Y68) and intein (Y69) in pSB1C3 (Y72+Y73) and pSB1K3 (Y70+Y71).</p> <p>The amount of digested insert used:</p> <p><u>Scaffold:</u> 20 fmol: 1,6 µl 50 fmol: 4 µl</p> <p><u>Intein</u></p>	iGEM2014_SOP015_v02_ligation

		<p>20 fmol: 1,8 µl 50 fmol 4,5 µl</p> <p><u>Backbones (10 fmol):</u> Y70: 1,6 µl Y71: 1,8 µl Y72: 1,33 µl Y73: 1,2 µl</p>	
15.08.2015	ADK	Transformation	iGEM2013_SOP0009_v0 1_TSB transformation
16.08.2015	TBA	<p>Colony-PCR on 3 colonies from each of the transformations. Colony-PCR showed appropriate bands for all except one, which was discarded.</p> <p>One successful colony of Intein and Scaffold in both backbones (4 in total) was selected for overnight culture.</p> <p>From Y70, 50 fmol colony 2 was selected for overnight culture.</p> <p>From Y71, 50 fmol colony 2 was selected for overnight culture.</p> <p>From Y72, 50 fmol colony 1 was selected for overnight culture.</p> <p>From Y73, 20 fmol colony 1 was selected for overnight culture.</p>	iGEM2013_SOP0021_v0 1_Colony PCR with MyTaq
17.08.2015	TBA, ADK	<p>Miniprep on the overnight Cultures and freezestock preparation. Miniprep stored as:</p> <ul style="list-style-type: none"> - R83: pSB1C3-initin - R84: pSB1C3-T18-linker-Scaffold-3xFLAG - R85: pSB1K3-initin 	iGEM 2014 SOP0019 - Miniprep

		- R86: pSB1K3-T18-linker-scaffold-3xFlag	
17.08.2015	JSP, TBA	<p>Received Target proteins (5 in total) from IDT along with Linker-Intein-Scaffold = LIS. Fast Digest of targets and LIS and appropriate backbones.</p> <p>Targets are linked to KT25, and were all digested with Xba1 + Pst1. Targets include; IL-2, Ybh, YbeY, Hfq and CsrA.</p> <p>LIS were digested with BamH1 + Pst1.</p> <p>Backbones were digestions of R32 (PcstA-RFP to fit targets) with Spe1 + Pst1 and R16 and R18 (to fit LIS) were digested with BamH1 + Pst1.</p> <p>Digestions were stored as:</p> <p>KT25_IL-2 = Y80 KT25_Ybh = Y81 KT25_YbeY = Y82 KT25_Hfq = Y83 KT25_CsrA = Y84 LIS = Y85 R32 (PcstA-RFP) = Y79 R16 (pSB1C3-UT18) = Y86 R18 (pSB1K3-UT18) = Y89</p>	iGEM2014_SOP0017_v01_ Fast digest
17.08.2015	ADK, TBA	<p>Ligations of the following:</p> <p>-LIS (Y85)/pSB1C3-UT18(Y86) -leuZ(Y63)/ pSB1C3-UT18(Y86) -leuZ(Y63)/pSB1C3-KT25(Y87) -leuZ(Y63)/pSB1K3-UT18(Y89) -leuZ(Y63)/pSB1K3-KT25 (Y88) -LIS(Y85)/pSB1K3-UT18(Y89)</p>	iGEM2014_SOP015_v02 _ligation

		<p>-T25_CsrA(Y84)/pSB1C3-PcstA(Y79)</p> <p>-T25_Hfq(Y38)/pSB1C3-PcstA(Y79)</p> <p>-T25_ybeY(Y82)/pSB1C3-PcstA(Y79)</p> <p>-T25_yhbJ(Y81)/pSB1C3-PcstA(Y79)</p> <p>-T25_iL2(Y80)/pSB1C3-PcstA(Y79)</p> <p>- Everything done according to the SOP</p>	
18.08.15	ADK	Transformation of ligations from yesterday, done according to the SOP	iGEM2013_SOP0009_v01_TSB transformation
24.08.15	JSP	Ligation of PcstA-RFP(Y79) and T25-IL2 (Y80) and T25-Hfq (Y83). 0, 20 and 50 fmol insert used. No alterations to SOP.	iGEM2014_SOP015_v02_ligation
25.08.15	JSP	Transformation of ligations from <i>24.08.2015</i> .	iGEM2013_SOP0009_v01_TSB transformation
26.08.15	JSP	Colony-PCR on transformations from 25.08.2015. No alteration to SOP.	iGEM2013_SOP0021_v01_Colony PCR with MyTaq
27.08.15	JSP	Miniprep on PcstA-RFP-T25-Hfq	Bio-Rad Quantum Prep: plasmid miniprep Kit
29.08.15	AC	Fast Digest of pSB1C3-T18-LIS. Cut with EcoI + PstI. Stored as Y105	iGEM2014_SOP0017_v01_Fast digest
	TBA	Ligation of Y105 into pSB1K3-backbone (Y104). Done	iGEM2014_SOP015_v02_ligation

		according to SOP and left overnight	
30-08.2015	EMT	Transformation done according to SOP	iGEM2013_SOP0009_v01_TSB transformation
31.08.2015	EMT	Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI + PstI. Ligation of digestion into pSB1K3-backbone done with 0, 20, and 50 fmol insert to 10 fmol backbone Transformation	iGEM2014_SOP0017_v01_Fast digest iGEM2014_SOP015_v02_ligation iGEM2013_SOP0009_v01_TSB transformation
01.09.2015	EMT	Ligation of digestion from 31.08.2015 into pSB1K3-backbone done according to SOP Transformation done according to SOP	iGEM2014_SOP015_v02_ligation iGEM2013_SOP0009_v01_TSB transformation
2.09.2015	EMT	Kolony PCR on pSB1K3-T18-LIS pSB1K3-LSx3FLAG done according to SOP	iGEM2013_SOP0021_v01_Colony PCR with MyTaq
3.9.2015	EMT	Miniprep of overnight Cultures done according to manufacturer's instructions	iGEM 2014 SOP0019 - Miniprep
5.08-2015	EMT	Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with XhoI. 100ng	iGEM2014_SOP0017_v01_Fast digest


		<p>done according to SOP</p> <p>Fast Digest of Library with XhoI</p> <p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library with a 1:10 scale (100fmol backbone to 1000fmol insert)</p> <p>Transformation done according to SOP</p>	<p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>
06.08.2015	EMT	<p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library (100fmol backbone to 1000fmol insert)</p> <p>Transformation: pSB1C3-T18-LSx3FLAG + Library pSB1C3-T18-LIS + Library done according to SOP</p>	<p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>
07.08.2015	EMT	<p>Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with XhoI.</p> <p>Fast Digest of Library 200ng with XhoI</p> <p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with Library</p> <p>Transformation pSB1C3-T18-LSx3FLAG + Library</p>	<p>iGEM2014_SOP0017_v01_Fast digest</p> <p>iGEM2014_SOP015_v02_ligation</p> <p>iGEM2013_SOP0009_v01_TSB transformation</p>

		pSB1C3-T18-LIS + Library done according to SOP	
08.09.2015	EMT	<p>Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI and PstI.</p> <p>Fast Digest of pSB1A2 with EcoI and PstI</p> <p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with pSB1A2 (0, 20 and 50fmol insert to 10fmol backbone)</p> <p>Transformation; pSB1A2 + T18-LIS pSB1A2 +T18-LSx3FLAG</p>	<p>iGEM2014_SOP0017_v01_ Fast digest</p> <p>iGEM2014_SOP015_v02 _ligation</p> <p>iGEM2013_SOP0009_v0 1_TSB transformation</p>
09.09.2015	EMT	<p>Ligation of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG with pSB1A2</p> <p>Transformation pSB1A2 + T18-LIS pSB1A2 +T18-LSx3FLAG</p>	<p>iGEM2014_SOP015_v02 _ligation</p> <p>iGEM2013_SOP0009_v0 1_TSB transformation</p>
10.09.2015	EMT	<p>Fast Digest of pSB1C3-T18-LIS and pSB1C3-T18-LSx3FLAG. Cut with EcoI and PstI.</p> <p>Fast Digest of pSB1A2 with EcoI and PstI</p> <p>Ligation of T18-LIS and T18-LSx3FLAG with pSB1A2</p> <p>Transformation; pSB1A2 + T18-LIS</p>	<p>iGEM2014_SOP0017_v01_ Fast digest</p> <p>iGEM2014_SOP015_v02 _ligation</p>


		pSB1A2 +T18-LSx3FLAG Koloni PCR on transformation from 09.09.2015	iGEM2013_SOP0009_v0 1_TSB transformation iGEM2013_SOP0021_v0 1_Colony PCR with MyTaq
11.09.2015	EMT	Koloni PCR on transformation from 10.09.2015	iGEM2013_SOP0021_v0 1_Colony PCR with MyTaq
12.09.2015	EMT, AC	Restriction Blot, done on pSB1A2-T18-LSx3FLAG, cuts made; -EcorI -EcorI+PstI -XhoI and uncut pSB1A2 Miniprep on T18-LSx3FLAG Fast Digest on pSB1A2-T18-LSx3FLAG with XhoI Ligation of digestion with Library Transformation of targets (T25-linker-CRSTA) into BTH101	iGEM 2014 SOP0019 - Miniprep iGEM2014_SOP0017_v01_ Fast digest iGEM2014_SOP015_v02_ligation iGEM2013_SOP0009_v0 1_TSB transformation

6. Results

Date (DD.MM.YY)	Picture	Comments
-----------------	---------	----------


11.08.2015		<p>R71 and R72 both contained BBa_J04450, which is the band visible at ~1000 bp.</p> <p>Bands just under 3000 was exercised from gel.</p>
12.08.2015		
13.08.2015		<p>Yesterdays transformations failed. might be due to mess up during plating out transformations</p>
14.08.15		
15.08.2015		

16.08.2015



Colony-PCR. Y71, 20 fmol colony 1 do not show appropriate band. This colony was discarded. Y71, 20 fmol colony 2 might have contained two colonies, one of which might have contained RFP instead of Intein. The second band corresponds to the length of RFP.

26.08.15



Only ligation of PcstA-RFP+T 25-Hfq was successful (#10). This was expected to be around 2700 nt. One of the bands showed a length of approx. 2700 nt. PcstA-RFP+T 25-IL2 was

		<p>expected be approx. 2600 nt, but were only approx. 1300 and 1600 nt.</p> <p>1-4: Y80-20 5-8: Y80-50 9-12: Y83-20 1-16: Y83-50</p>
30.08.2015		Transformation Failed
02.09.2015		Transformation of T18-LIS and T18-LSx3FLAG succeeded. Checked with colony PCR
10.09.2015		Colony PCR was inconclusive
11.09.2015		Colony PCR showed that ligation of pSB1A2-T18-LSx3FLAG was successful, and pSB1A2-T18-LIS failed

7. Appendices

