

# iGEM 2015 – Microbiology – BMB – SDU

**Project type:** Interlab study.

**Project title:** iGEM interlab study of GFP.

**Sub project:** Measurement of GFP

**Creation date:** 23/6 2015

**Written by:** CEM, AC

**Performed by:** AC

## 1. SOPs in use:

iGEM 2014 SOP0009 - TSB Transformation.

iGEM 2013 SOP0021 - Colony PCR with MyTag

iGEM 2013 SOP0005 - Freeze Stock

iGEM 2014 SOP0019 - Miniprep

iGEM 2014 SOP0017 - Fast Digest

iGEM 2014 SOP0015 - Ligation

## 2. Purpose:

**Measurement of GFP under influence of different promoters.**

## 3. Overview

Date (DD-MM-YY)	Person(s) (initials)	Experiments	SOPs
23-06-15	AC	Transformation	iGEM 2014 SOP0009
24-06-15	AC	Colony PCR	iGEM 2013 SOP0021
25-06-15	AC	Freeze stock	iGEM 2013 SOP0005
25-06-15	AC	Miniprep	iGEM 2014 SOP0019
01-07-15	AC	Fast digest	iGEM 2014 SOP0017
01-07-15	AC	Ligation	iGEM 2014 SOP0015
02-07-15	AC	Transformation	iGEM 2014 SOP0009
03-07-15	AC	Transformation	iGEM 2014 SOP0009
04-07-15	AC	Freezestock	iGEM 2013 SOP0005

<b>04-08-15</b>	<b>AC</b>	<b>Colony PCR</b>	<b>iGEM 2013 SOP0021</b>
<b>Start from scratch</b>			
<b>05-08-15</b>	<b>AC</b>	<b>Transformation</b>	<b>iGEM 2014 SOP0009</b>
<b>06-08-15</b>	<b>AC</b>	<b>Miniprep</b>	<b>iGEM 2014 SOP0019</b>
<b>06-08-15</b>	<b>AC</b>	<b>Fast Digest</b>	<b>iGEM 2014 SOP0017</b>
<b>06-08-15</b>	<b>AC</b>	<b>Ligation</b>	<b>iGEM 2014 SOP0015</b>
<b>07-08-15</b>	<b>AC</b>	<b>Transformation</b>	<b>iGEM 2014 SOP0009</b>
<b>09-08-15</b>	<b>AC</b>	<b>Colony PCR</b>	<b>iGEM 2013 SOP0021</b>
<b>10-08-15</b>	<b>AC</b>	<b>Freeze stock</b>	<b>iGEM 2013 SOP0005</b>
<b>10-08-15</b>	<b>AC</b>	<b>Miniprep</b>	<b>iGEM 2014 SOP0019</b>

#### **4. Materials required.**

##### **Materials in use**

<b>Name</b>	<b>Components (Concentrations)</b>	<b>Manufacturer / Cat. #</b>	<b>Room</b>	<b>Safety considerations</b>
<b>EcoRI</b>		<b>Roche: Fast Digest</b>		
<b>XbaI</b>		<b>Roche: Fast Digest</b>		
<b>SpeI</b>		<b>Roche: Fast Digest</b>		
<b>Pst</b>		<b>Roche: Fast Digest</b>		
<b>FACS</b>				
<b>PBS</b>				
<b>LB medium</b>				
<b>LB agar medium</b>				

<b>PCR machine</b>				
<b>Equipment for gel-electrophoresis</b>				




## 5. Experiment history

<b>Date (DD-MM-YY)</b>	<b>Person(s) (initials)</b>	<b>Alterations to SOPs and remarks to experiments</b>	<b>SOPs</b>
<b>23-06-15</b>	<b>AC</b>	Plasmids from plate 1 well 20K, 22K and 22A transformed into MG1655. Plasmid from plate 4 well 21J	iGEM 2014 SOP0009
<b>24-06-15</b>	<b>AC</b>	Colony PCR, of 2 cultures per plate	iGEM 2013 SOP0021
<b>25-06-15</b>	<b>AC</b>	Freeze stock	iGEM 2013 SOP0005
<b>25-06-15</b>	<b>AC</b>	Miniprep	
<b>01-07-15</b>	<b>AC</b>	Fast digest, GFP cut with EcoI and Pst. Rest cut with SpeI and Pst.	iGEM 2014 SOP0017
<b>01-07-15</b>	<b>AC</b>	Ligation with 0 fmol insert, and 50 fmol insert.	iGEM 2014 SOP0015
<b>02-07-15</b>	<b>AC</b>	Transformation into Top10	iGEM 2014 SOP0009
<b>04-07-15</b>	<b>AC</b>	Colony PCR (No gel)	iGEM 2013 SOP0021
<b>05-07-15</b>	<b>CEM</b>	Gel on the colony PCR from 04-07-15	
<b>06-07-15</b>	<b>AC</b>	Miniprep	iGEM 2014 SOP0019
<b>06-07-15</b>	<b>AC</b>	Transformation into MG1655	iGEM 2014 SOP0009
<b>08-07-15</b>	<b>AC</b>	Transformation into MG1655	iGEM 2014 SOP0009

<b>09-07-15</b>	<b>AC</b>	Transformation into MG1655	iGEM 2014 SOP0009
<b>29-07-15</b>	<b>AC</b>	Transformation into MG1655	iGEM 2014 SOP0009
<b>30-07-15</b>	<b>AC</b>	Miniprep of GFP and plasmids containing different promoters	iGEM 2014 SOP0019
<b>30-07-15</b>	<b>AC</b>	Fast Digest of GFP and promoters	iGEM 2014 SOP0017
<b>30-07-15</b>	<b>AC</b>	Ligation of plasmids with GFP as insert	iGEM 2014 SOP0015
<b>31-07-15</b>	<b>AC</b>	Transformation into Top10	iGEM 2014 SOP0009
<b>01-08-15</b>	<b>AC</b>	Colony PCR	iGEM 2013 SOP0021
<b>02-08-15</b>	<b>AC</b>	Colony PCR	iGEM 2013 SOP0021
<b>06-08-15</b>	<b>AC</b>	Miniprep	iGEM 2014 SOP0019
<b>13-08-15</b>	<b>AC</b>	FACS done according to iGEM protocol and worksheet	
<b>19-08-15</b>	<b>AC</b>	FACS done according to iGEM protocol and worksheet	
<b>25-08-15</b>	<b>AC</b>	FACS done according to iGEM protocol and worksheet	

## 6. Results

<b>Date (DD-MM-YY)</b>	<b>Picture</b>	<b>Comments</b>
<b>03-07-15</b>		Transformation inconclusive at this point.

<p><b>05-07-15</b></p>		<p>The Gele showed that the plasmid was around 1200 bp long.</p>
<p><b>07-07-15</b></p>		<p>Transformation unsuccessful (Wrong plates)</p>
<p><b>09-07-15</b></p>		<p>Transformation unsuccessful (Wrong plates)</p>
<p><b>01-08-215</b></p>		<p>Colony PCR confirmed Ligation for 2 out of three plasmids</p>
<p><b>02-08</b></p>		<p>Colony PCR confirms last plasmid</p>
<p><b>06-08-2015</b></p>		<p>Miniprep made from each promoter</p>
		<p>Sequencing</p>
<p><b>25-08-2015</b></p>		<p>FACS Measurement</p>
<p><b>26-08-2015</b></p>		<p>Two constructs confirmed through sequencing. Construct 1 (J23101+I13504) showed to have a faulty sequence, but did show fluorescence when analyzed through FACS</p>

## **7. Appendices**