

Ligation Protocol (Promega)

Materials:

- Vector DNA (~100ng)
- Insert DNA
- Ligase 10x buffer
- T4 DNA ligase
- Nuclease-free water

Protocol:

Promega recommends using a 1:1, 1:3 or 3:1 molar ratio of vector:insert DNA when cloning a fragment into a plasmid vector. The ng of insert can be obtained by using the below equation.

$$\frac{\text{vector (ng)} \times \text{insert (kb)}}{\text{vector (kb)}} \times \text{molar ratio of } \frac{\text{insert}}{\text{vector}} = \text{insert (ng)}$$

1. Assemble the following reaction in a sterile microcentrifuge tube:

Vector DNA	100 ng
Insert DNA	From equation
Ligase 10x Buffer	1 µl
T4 DNA Ligase (Weiss units)	0.1-1
Nuclease-Free Water to final volume of	10 µl

1. Incubate the reaction at room temperature for 3 hours, or at 4°C overnight.
2. Transform with 2 µl of product.