

Protocol 10: RT-PCR

1. Material

- RNA
- RT Buffer
- dNTPs (2.5mM each)
- dNTPs (10 mM each)
- oligo dT 18mer (100 μ M)
- DEPC-H₂O
- RevertAid H Minus RT, Fermentas
- RNasin (Ribolock RI, Fermentas)
- Taq-Pol
- Primers (10 μ M)
- PCR Tubes

2. Instruments

- Thermocycler

3. Experimental procedure

-Use a PCR tube for sample preparation.

-Prepare a master mix

Final volume 25 μ l

X μ l	1 μ g RNA
5 μ l	5xRT Buffer
2.5 μ l	dNTPs (2.5mM each)
0.5 μ l	oligo dT 18mer (100 μ M)
Add to 23 μ l	DEPC-H ₂ O (bidest, sterile)

Set up program:

-incubate 5min at 70°C (denaturing)

-incubate for 5min at 37°C (primer annealing)

-Pause

-Add 1 μ l RevertAid H Minus RT, Fermentas 200U/ μ l and 1 μ l RNasin (Ribolock RI, Fermentas 40 U/ μ l)

-incubate 1 hour at 42°C

-incubate 15 min at 70°C (inactivating the enzymes, thus stopping the reaction)

-store cDNA at -20°C or start directly a PCR

PCR reaction

Sample preparation on ice

cDNA	1µl
10xPCR buffer	5µl
Primer1 (10µM)	0.5µl
Primer2 (10µM)	0.5µl
dNTPs (10mM each)	1µl
Taq-Pol (xU/µl)	0.25µl
H2O	41µl

Set up your PCR program, eg:

1. 94°C 5min
2. 94°C 30sec
3. 58°C 30sec (temperature depending on your particular primer pair)
4. 72°C 90sec, repeat step 2 to 4 30 times
5. 72°C 10min (extension)
6. 4°C pause