

Determination of CCH2 rate:

Materials Needed:

250 mL Erlenmeyer flasks
M9 minimal salt media
Chloramphenicol (25mg/ml)
1M IPTG
200 mM 6-CNA dissolved in methanol
250 mL Erlenmeyer flasks with 50 mL LB
50 mL Falcon Tubes
Spectrophotometer
Cuvette
1.5mL Eppendorf tubes
Centrifuge
GC/MS
Pipet

M9 minimal salt media components:

33.9g/L Na₂HPO₄
15g/L KH₂PO₄
5g/L NH₄Cl
2.5g/L NaCl

Protocol:

- 1) A starter culture is prepared from an isolated colony of DH5- α E.coli transformed with a CCH2 vector in 50mL of LB with 25ug/mL of Chloramphenicol overnight on the shaker at 37°C.
- 2) A 'no insert' colony is used as a control.
- 3) 500 μ L of the starter culture is used to inoculate 50mL of LB with 25ug/mL of chloramphenicol in a 250 mL Erlenmeyer flask and incubated on the shaker at 37°C.
- 4) Cultures are grown until OD₆₀₀= 0.6-0.8, then put on ice for 20 minutes and stirred every couple of minutes.
- 5) After being cooled, 50 μ L of 1M IPTG is added and the Erlenmeyer stirred.
- 6) The cultures are grown overnight on the shaker at 25°C.
- 7) Cells are pelleted in 50 mL falcon tubes at 4°C at 4000 rpm for 20 min.
- 8) The supernatant is discarded and the cultures are resuspended in 40 mL of M9 minimal salt media to wash the cells
- 9) Cells are pelleted again, as above, and the supernatant is discarded
- 10) The cells are resuspended in 25mL of M9 minimal salt media.
- 11) The duplicate samples are combined into one 50mL falcon tube to normalize the amount of cells in each fraction.
- 12) 25mL of 50mL from the fractions is transferred into two separate, empty 250 mL Erlenmeyer flasks.
- 13) 20 μ L of 200mM 6-CNA is added to all samples and were grown on the shaker for four hours at 30°C.

- 14) An initial sample is collected at time zero. Subsequently samples are collected every 15 minutes over a period of two hours. Thereafter, samples are collected every 30 minutes for another two hours.
- 15) Collected samples were centrifuged for 2.5 minutes at 15000rpm at 4°C. The supernatants of the samples were transferred into 1.5mL Eppendorf tubes and stored at -80°C .
- 16) The samples are extracted via a liquid-liquid extraction to isolate 6-CNA and 6-HNA and ran on the GC-MS [\[1,2\]](#).

Other Protocols:

[\[1\]](#) Gas Chromatography-Mass spectroscopy conditions to run 6-Chloronicitinoic Acid and 6-hydroxynicitinoic Acid and sample preparation

[\[2\]](#) Liquid-Liquid Extraction for isolation of 6-Chloronicitinoic Acid and 6-Hydroxynicitinoic Acid