



**UNSW**  
AUSTRALIA

BABS UNSW iGEM Lab Protocol



Procedure	Name		Restriction Enzyme Digestion			
	Description		Enzyme Digestion of DNA			
Document	Name	Isabelle Capell-Hattam	Date	2/07/15	Version	1
Requirements	Time		1.5 - 2 hours			
	PPE		Gloves, Labcoat			
	Equipment		Pipettes and Tips Heat Block 1.5 mL tubes Ice bucket			
	Materials		NEB Buffer 2 BSA EcoR1 Pst1 Spe1 Xba1 MiliQ H2O Dpn1 (optional)			
Digestion of a single fragment for insertion into backbone						
Step 1	Add into each tube: <ul style="list-style-type: none"> <li>● 250 ng DNA</li> <li>● 2.5 <math>\mu</math>L of NEB Buffer 2</li> <li>● 0.5 <math>\mu</math>L BSA</li> <li>● 0.5 <math>\mu</math>L EcoR1</li> <li>● 0.5 <math>\mu</math>L Pst1</li> </ul> Make the total volume in the tube up to 20 $\mu$ L with MiliQ H2O					
Step 2	Incubate at 37°C for 30 minutes					
Step 3	Heat inactivate at 80°C for 20 minutes					
Step 4 (optional)	Run digested products on a gel to verify size					
Step 5 (optional)	Ligate fragments together					

Digestion of 2 fragments for directional assembly into a backbone	
Step 1	Add into each tube: <ul style="list-style-type: none"> <li>● 250 ng DNA</li> <li>● 2.5 <math>\mu</math>L of NEB Buffer 2</li> <li>● 0.5 <math>\mu</math>L BSA</li> </ul> Make the total volume in the tube up to 19 $\mu$ L with MiliQ H <sub>2</sub> O
Step 2	In the backbone tube add: <ul style="list-style-type: none"> <li>● 0.5 <math>\mu</math>L EcoR1</li> <li>● 0.5 <math>\mu</math>L Pst1</li> </ul> In the fragment 1 tube add: <ul style="list-style-type: none"> <li>● 0.5 <math>\mu</math>L EcoR1</li> <li>● 0.5 <math>\mu</math>L Spe1</li> </ul> In the fragment 2 tubes add: <ul style="list-style-type: none"> <li>● 0.5 <math>\mu</math>L Xba1</li> <li>● 0.5 <math>\mu</math>L Pst1</li> </ul>
Step 3	Incubate at 37°C for 30 minutes
Step 4	Heat inactivate at 80°C for 20 minutes
Step 5 (optional)	Run digested products on a gel to verify size
Step 6 (optional)	Ligate fragments together
Notes	This protocol was adapted from the 2014 iGEM HQ protocol. Enzymes should be kept on ice at all times 0.5 $\mu$ L of Dpn1 can be added to the linearized backbone digest to prevent background colonies from being observed after ligation and transformation
Version History	