

## **Week 4: 29. June 2015- 3. July 2015**

### **29. June 2015**

#### 1) Gel-Extraction of *mmoD* and pET-30

- Use the samples frozen at 26. June 2015
- Use the Wizard® SV Gel and PCR Clean-Up System Kit (Promega) and follow the provided manual

#### 2) Ligation of *mmoD* into pET-30

- Pipetting scheme according to the Protocol Ligation with T4 DNA Ligase (NEB)
- Control Reaction: pET-30 Selfligation
- Incubate the reaction 30 minutes at room temperature

#### 3) Transformation of *E. coli* TOP10 cells with ligation reaction pET-30+*mmoD* and pET-30 selfligated

- Add 5 µl of Ligation reaction to chemically competent *E. coli* TOP10 cells.
- Following the Protocol for Transformation of chemical competent *E. coli* cells
- Plate on LB+Kan [50 µg/ml] and incubate overnight at 37 °C

#### 4) Isolation of genomic DNA from *Methylococcus capsulatus*


- Use NucleoSpin® Tissue Kit (Machery-Nagel) and follow the provided manual

### **30. June 2015**

- The Ligation does seem to work, the control plate, containing the transformed *E. coli* TOP10 with the selfligated pET-30 had an equal amount of colonies compared to the transformed *E. coli* TOP10 with the pET-30+*mmoD* Ligation

#### 1) Digest of *mmoB*, *mmoX* and pET-30 with NdeI/EcoRI

- Pipetting scheme according to the protocol restriction digest
- First adding NdeI and incubate the Reaction at 37 °C for 1 hour, afterwards heat inactivation of NdeI at 65 °C for 2 minutes. Adding 1 µl of EcoRI and incubate the digest again for 1 hour at 37 °C
- Verify 10 µl of the digestion on agarose gel
- Control: undigested pET-30



**Figure 1: Digest of *mmoX*, *mmoB* and pET-30 with NdeI/ EcoRI.** 10  $\mu$ l of Digest were checked on 1% (w/v) agarose gel. Expected sizes: pET-30-5269 bp, *mmoX*- 1605 bp, *mmoB*- 447 bp. Undigested pET-30 was included as negative control (NC). As ladder (L) 1kB Ladder (NEB) was used.

## 2) Purification of digested *mmoB*, *mmoX* and pET-30


- Using Wizard® SV Gel and PCR Clean-Up System (Promega) and following the provided manual

## 3) Ligation of *mmoX* and *mmoB* into pET-30

- Pipetting scheme according to the Protocol Ligation with T4 DNA Ligase (NEB)
- Control Reaction: pET-30 Selfligation
- Incubate the reaction at 16 °C overnight

## 4) Amplification of *mmoG* adding Restriction sites (NdeI/EcoRI) for cloning into pET-30

- Pipetting scheme and PCR program according to PCR with Phusion-HF DNA Polymerase protocol
- Primer: *mmoG*\_E1/E2,
- Template: PCR Product *mmoG* for TOPO cloning (26. June 2015)



**Figure 2: Amplification of *mmoG* adding restriction sites for cloning into pET-30.** 20  $\mu$ l of PCR were checked on 1% (w/v) agarose gel. Expected size: *mmoG*. As ladder (L) 1kB Ladder (NEB) was used.

## 1. July 2015

### 1) Transformation of *E. coli* TOP10 cells with pET-30+mmoX, pET-30+mmoB and pET-30 selfligated


- Add 5 µl of ligation reaction to chemically competent *E. coli* TOP10 cells.
- Following the protocol for transformation of chemical competent *E. coli* cells
- Plate on LB+Kan [50 µg/ml] and incubate overnight at 37 °C

## 2. July 2015

- The ligation was successful, control reaction with selfligated pET-30 showed less clones compared to the ligation reactions

### 1) Colony-PCR to screen for clones containing the pET-30+mmoX or pET-30+mmoB

- Constructs: pET-30+mmoX, pET-30+mmoB
- Pipetting scheme and PCR program according to PCR with Taq-DNA Polymerase Protocol
- Primer: mmoB\_E1/E2, mmoD\_E1/E2
- Check 6 Clones per Construct
- Positive Control: add 1 µl of pSC1B3+mmoX or pSC1B3+mmoB
- Negative Control: add 1 µl MilliQ Water



**Figure 3: Colony-PCR to identify clones containing pET-30+mmoX or pET-30+mmoB expression constructs.** Numbers 1-6 determine the checked clone. As positive control (PC) 1 µl pSC1B3+mmoX or pSC1B3+mmoB was added. As negative control (NC) 1 µl MilliQ Water was added. 10 µl of PCR were analyzed on 1 % (w/v) agarose gel. Expected sizes: mmoB-453 bp, mmoX-1611 bp. As ladder (L) 1 kb Ladder (NEB) was used.

- The pET-30+mmoB expression plasmid seems to be successful cloned

1) Inoculation of liquid culture for plasmid isolation of pET-30+*mmoB* and pET-30 selfligated

- Inoculate clone 2 and 4 with 5 ml LB+Kan [50 µg/ml] and incubate overnight at 37 °C shaking at 220 rpm.
- Inoculate one of the clones of the pET-30 selfligation with 5 ml LB+Kan [50 µg/ml] and incubate overnight at 37 °C shaking at 220 rpm.

**3. July 2015**

1) Plasmid Isolation of pET-30+*mmoB*

- Use the QIAprep Spin Miniprep Kit (Qiagen) and follow the provided manual.