## Test of detection system in Saccharomyces cerevisiae

## Extraction of h+ pheromones from *S. pompe*:

*S. pombe* L972 +h was grown in YPD media at 30  $^{\circ}$ C for 3 weeks. The pheromones was extracted pelleting the cells and filtrating the supernatant through filter with 0.2  $\mu$ M pore size. To increase the concentration of the h+ phermones the solution was run in a concentration with cold trap for 5,5 h. The concentration was increased by approximately 20 times.

## Test of detection system:

- 1. A colony of *S. cerevisiae* CEN.PK2 with integrated construct 4 was inoculated in 5 ml YPD media overnight. Wild type CEN.PK2 and a strain expressing RFP was used as negative and positive control.
- 2. The OD was measured after overnight preculture and the cells was centrifuged at 1100 rcf for 5 min.
- 3. The pellet was dissolved in 5 ml YPD media and transferred to a shake flask. The solution was diluted with YPD to a OD of 0,4. The cells was inoculated at 30 °C for 2h.
- 4. 0, 50, 115 and 230  $\mu$ L concentrated pheromone was added to 1 ml cell suspension of the negative control and the colony with C4. No pheromones was added to the positive control. The cell was once again incubated at 30 °C for 2h to allow expression of RFP.
- 5. The cells were pelleted and washed with once with  $dH_20$ . The pellet was dissolved in  $70\mu l$   $dH_20$  and 3  $\mu L$  of the solution was used for studying of the detection system with a fluorescence microscope.
- 6. Used 1 second of exposure time for RFP measurement. Used GFP exposure as a measure of inviable cells.