

Chemocompetent *Escherichia coli* cells

This protocol is from the iGEM team Aix-Marseille 2014: http://2014.igem.org/Team:Aix-Marseille/Protocol:chemocompetent_ecoli_cells

Chemoconpetent *Escherichia coli* cells

1. Streak *E.coli* cells on an LB plate with or without antibiotics.
2. Allow cells to grow at 37°C overnight.
3. Place three or four colonies in 5 mL LB media (+ antibiotic selection if necessary), grow overnight at 37°C.
4. Dilute the preculture to OD600 = 0,05 in 100 mL LB.
5. Allow cells to grow at 37°C (250 rpm), until OD600= 0.4 (about 2-3 hours).
6. Transfer cells to 2 centrifuge bottles (50 mL), and place cells on ice for 20 minutes.
7. Centrifuge cells in Sorval GSA rotor at 4°C for 10 minutes at 5 000 rpm.
8. Subsequent resuspensions may be done in the same bottle. Cells must remain cold for the rest of the procedure: Transport tubes on ice and resuspend on ice in the cold room. Discard the supernatant and resuspend cells in 1/2 Vol of cold 50 mM CaCl₂. Incubate on ice for 20 minutes.
9. Centrifuge cells using Sorval RT6000B rotor at 4°C for 10 minutes at 5 000 rpm.
10. Discard the supernatant and resuspend cells (by pipetting) in 1/20 Vol cold 50 mM CaCl₂ containing 15% glycerol. Transfer 500 µL of cells into (1.5 mL) Ependorff tubes placed on ice. Cells stored at -80°C can be used for transformation for up to ~6 months.

Note: Through the process, cells should be treated with care. No vortexing or excess pipetting should be performed, especially when the cells have been resuspended in CaCl₂ because lysis will result, decreasing the amount of competent cells. Also, depending on the density of the cells, higher or lower volumes CaCl₂ can be used to increase the concentration of cells per tube.

Transformation

11. Add 2 µL of 0,5, 10 and 50 ng/L of Transformation Efficiency Kit DNA to 100 µL of competent cells.
12. Incubate the mixture on ice for 35 minutes.
13. Transfer the reaction to a 42°C water bath for exactly 2 minutes.
14. Incubate on ice for 5 minutes.

15. Add 900 μ L of LB medium to each tube and incubate at 37°C for 1 hour to allow cells to recover and express the antibiotic resistance marker.
16. Spread the appropriate quantity of cells (50 to 200 μ l) on selective media. Store the remaining cells at 4°C.
17. Incubate all plates overnight at 37°C (agar side up).