

2015SZU-iGEM PROJRCT PROTOCOL

一、ACKRS (1260bp)

ACKRS needs to be PCR amplified with the primer,so that it could have the required enzyme cutting sites.

1、PCR clone

1.1、PCR try

1.1.1、PCR trial system 1:

ACKRS template :	1 μ L
ACKRS (王皓 EcoR1) :	0.25 μ L
ACKRS (王皓 Pst1) :	0.25 μ L
H2O:	3.5 μ L
Mix:	5 μ L
<hr/>	
total:	10 μ L

Set three temperature gradients:

- ① 95°C: 5min → 94°C: 30s → 48°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞
- ② 95°C: 5min → 94°C: 30s → 55°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞
- ③ 95°C: 5min → 94°C: 30s → 58°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞

After the PCR, run agarose gel electrophoresis.If the bright band was in the 1260bp position, indicating that PCR was successful, you can directly enter the step 1.2.1. If not you should try 1.1.2 system.

1.1.2、PCR trial system 2:

ACKRS template :	1 μ L
ACKRS (王永益 EcoR1) :	0.25 μ L
ACKRS (王永益 Pst1) :	0.25 μ L
H2O:	3.5 μ L
Mix:	5 μ L
<hr/>	
total:	10 μ L

Set three temperature gradient:

- ④ 95°C: 5min → 94°C: 30s → 48°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞
- ⑤ 95°C: 5min → 94°C: 30s → 55°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞
- ⑥ 95°C: 5min → 94°C: 30s → 58°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞

After PCR should run the agarose gel electrophoresis.If the light band was in the 1260bp position, indicating that PCR was successful,used 1.1.3 system to PCR again.

1.1.3、PCR trial system 3:

PCR trial system 1 product:	1 μ L
ACKRS (王皓 EcoR1) :	0.25 μ L
ACKRS (王皓 Pst1) :	0.25 μ L
H2O:	3.5 μ L
Mix:	5 μ L
<hr/>	
total:	10 μ L

Set three temperature gradient:

- ⑦ 95°C: 5min → 94°C: 30s → 48°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞
- ⑧ 95°C: 5min → 94°C: 30s → 55°C: 30s → 72°C: 2min → 72°C: 7min → 4°C: ∞

2015SZU-iGEM PROJRCT PROTOCOL

- ⑨ 95°C: 5min→ 94°C: 30s→ 58°C: 30s→ 72°C: 2min→ 72°C: 7min→ 4°C: ∞

After Pcr,should run agarose gel electrophoresis.If the light band was in the 1260bp position, indicating that PCR was successful, you can directly enter the step 1.2.2.

1.2、PCR amplification

1.2.1 PCR amplification system 1:

ACKRS plasmid template:	1 μ L
ACKRS (王皓 EcoR1) :	0.5 μ L
ACKRS (王皓 Pst1) :	0.5 μ L
H2O:	13 μ L
Mix:	15 μ L
<hr/>	
total:	30 μ L

- ① The system was formulated with 5 tubes, that was, 30ul*5.
- ② Temperature setting should be and the same with the temperature in step 1.1.1 .
- ③ After PCR, if the bright band was in the position of 1260bp,indicated the PCR was successful, and next should purify the PCR product.
- ④ Made marks: PCR ACKRS purification + date, and entered step 2.

1.2.2 PCR amplification system 2:

Taking 1.1.2 step PCR product as a template, made a large PCR amplification.

PCR trial system 1 product:	1 μ L
ACKRS (王皓 EcoR1) :	0.5 μ L
ACKRS (王皓 Pst1) :	0.5 μ L
H2O:	13 μ L
Mix:	15 μ L
<hr/>	
total:	30 μ L

- ① After PCR, if the light band was in the position of 1260bp,indicated the PCR was successful, and the next step should purify the PCR product.
- ② Made marks: PCR ACKRS purification + date, and enter step 2.

2、Enzyme digestion

ACKRS PCR purification product:	3 μ L
Xba1:	1 μ L
Pst1:	2 μ L
Tango:	1 μ L
H2O:	3 μ L
<hr/>	
total:	10 μ L

- ① The enzyme digestion system was formulated with 10 tubes, which were digested with 10 L*10.
- ② Enzyme digestion condition: 37 C, 3h above (digestion overnight is best , also can digest all day.It is good for ligation overnight by T4 ligase)
- ③ After enzyme digestion, the system was merged 10 tubes into one tube, and purified after running agarose gel electrophoresis,which showed the bright band in 1260bp.
- ④ Made marks:ACKRS (XP) digestion and purification+ date, and than entered step 3.

3、Ligation

Ligate ACKRS with hUPII and pSB1C3

ACKRS (XP) :	3 μ L
hUPII (SP) :	3 μ L
T4 ligase :	1 μ L
10*T4 buffer :	2 μ L
H2O :	1 μ L
<hr/>	
Total:	10 μ L

- ① The ligation system needed to be formulated with 5 tubes, and each tube was 10 ul.
- ② Ligation condition: 16, 3h (best ligation overnight).
- ③ After the ligation, could directly enter step 4.

4、Transformation

- ① Packed a box of ice. Took a tube of competent cells(100ul) from -80 degrees refrigerator and placed in the ice to melt (can not be placed at room temperature, nor a violent shock, otherwise the competent cells would be dead).
 - ② When competent cells melted completely(about 10min), added all the ligation product (10 ul) the competent cells, and gently blew. (not a violent shock, can not come out of the ice box, it is necessary to carry on the ice.)
 - ③ Placed 30min on ice.
 - ④ Set the competent cells(110ul) in sponge and placed in 42 $^{\circ}$ C water to bath 90s (not too long), immediately removed and placed in the ice, ice bath for five minutes.
 - ⑤ Lighted alcohol burner in the super clean bench, added 890ul LB to competent cells (the super clean bench before using should be disinfected by ultraviolet,the LB before using and covering are required to the long time firing.)
 - ⑥ Put the competent cells which added LB on the table concentrator (set in the sponge) ,150rpm shaking an hour(table password is: 3).
 - ⑦ An hour later,took out the competent cells,and centrifugated with 12000rpm 3min.
 - ⑧ Took a chloromycetin(Chl) LB Plating. Removed 900ul LB of competent cells,which was centrifugated, leaving 100ul competent cells,and blew 100ul competent cells with transfer liquid gun,and than added it to the Chl LB Plating. Used glass spreader (before using,needed to be cool after alcohol burner firing,about 5min) to flat the bacterial liquid gently on the LB Plating. Sealed the sealing film, and finally made good marks: hUPII+AckRS+ date, placed in a 37 $^{\circ}$ C constant temperature incubator for 2 days .
- Two days later,entered step 5.

5、Picking bacterial colonies and Culture

- ① Took ten 1.5ml EP tubes from the lunch box in the super clean bench, numbered 1-10, each EP tube added in 1ml LB liquid medium and 1ul Chl. Picked 10 bacterial colonies and respectively put into the 10 EP tubes.

- ② Set the 10 EP tubes into the sponge, put it in the table concentrator, shaking with 220rpm, until the medium became turbid.
- ③ Entered step 6.

6、PCR verification

- ① Put the 10 EP tubes of step 5 in the centrifuge to centrifugate in 12000rpm 3min.
- ② Preparation of the following PCR system:

ACKRS (up)	:	0.25 μ L
ACKRS (down)	:	0.25 μ L
hUPII (up)	:	0.25 μ L
hUPII (down)	:	0.25 μ L
H2O:		4 μ L
Mix :		5 μ L
<hr/>		
total:		10 μ L

- ③ The PCR system needed to be formulated with 10 tubes, and numbered them 1-10.
- ④ Sucked up about 1ul bacterial precipitation to add to the corresponding PCR tubes with the small transfer liquid gun in the super clean bench, and PCR at the optimum temperature of step 1.1.
- ⑤ Took 5ul from each of 10 tubes and verified by running agarose electrophoresis. If there were bright bands at 1260bp (AckRS) and 355bp (hUPII), needed to transfer the bacteria to a Erlenmeyer flask to take a large amount of culture (adding about 20ml LB, and the addition 20ul of chloromycetin (Chl)).
- ⑥ Extracted plasmid, and verified by single and double enzyme digestion.
- ⑦ Prepared the following enzyme digestion system:

Single enzyme digestion:		double enzyme digestion:	
plasmid:	5 μ l	plasmid:	5 μ l
Pst1:	1 μ l	EcoR1:	1 μ l
10*Buffer O:	2 μ l	Pst1:	1 μ l
H2O:	12 μ l	10*Buffer O:	2 μ l
		H2O:	11 μ l
<hr/>		<hr/>	
total:	20 μ l	total:	20 μ l


- ⑧ If the single and double digestion could cut out the corresponding size of the bright band, it indicated the BioBrick was successfully constructed.

二、tRNA (CUA) (72bp)

The full-length of the gene of tRNA is 72bp, we got it by synthesis, and its left end has EcoRI and XbaI sites, the right end contains spel and pstI sites.

1、Enzyme digestion

tRNA:	:	5 μ L
EcoR1:	:	1 μ L
Pst1:	:	1 μ L
10*Buffer O :		2 μ L
H2O:		11 μ L
<hr/>		
total:		20 μ L



2015SZU-iGEM PROJRCT PROTOCOL

- ① The digestion system needed to be formulated with 2 tubes, that is 20ul*2.
- ② Enzyme cutting conditions: 37 C, 3h above (digested overnight is the best, also can digest all day.)
- ③ Purified after the digestion.Enter step 2.

2、Ligation

Ligated tRNA to pSB1C3

tRNA (EP) :	5 μ L
PSB1C3 (EP) :	1 μ L
T4 ligase :	1 μ L
10*T4 buffer:	2 μ L
ddH2O:	1 μ L
<hr/>	
total:	10 μ L

- ① The connecting system needed to be formulated with 5 tubes, that is 10ul*5., which are connected together.
- ② Connection condition: 16, 3h (best ligated overnight).
- ③ After the ligation,did not need to run agarose electrophoresis, and entered step 3.

3、Transformation

- ① Packed a box of ice. Took a tube of competent cells(100ul) from -80 degrees refrigerator and placed in the ice to melt (can not be placed at room temperature, nor a violent shock, otherwise the competent cells would be dead).
- ② When competent cells melted completely(about 10min), added all the ligation product (10 ul) the competent cells, and gently blew. (not a violent shock, can not come out of the ice box, it is necessary to carry on the ice.)
- ③ Placed 30min on ice.
- ④ Set the competent cells(110ul) in sponge and placed in 42℃ water to bath 90s (not too long), immediately removed and placed in the ice, ice bath for five minutes.
- ⑤ Lighted alcohol burner in the super clean bench, added 890ul LB to competent cells (the super clean bench before using should be disinfected by ultraviolet,the LB before using and covering are required to the long time firing.)
- ⑥ Put the competent cells which added LB on the table concentrator (set in the sponge) ,150rpm shaking an hour(table password is: 3).
- ⑦ An hour later,took out the competent cells,and centrifugated with 12000rpm 3min.
- ⑧ Took a LB plating withchloromycetin(Chl).Removed 900ul LB of competent cells,which was centrifugated,leaving 100ul competent cells,and blew 100ul competent cells by transfer liquid gun,and than added it to the Chl LB Plating. Used glass spreader (before using,needed to be cool after alcohol burner firing,about 5min) to flat the bacterial liquid gently on the LB Plating. Sealed the sealing film, and finally made good marks: hUPII+AckRS+ date, placed in a 37℃ constant temperature incubator for 2 days .
- ⑨ Two days later,entered step 5 and step 6 of “一、ACKRS (1260bp)”