

## **Week 12: 24. August 2015 – 28. August 2015**

### **24. August 2015**


#### 1) Expression of MEDH2 into *E. coli* BL21 for expression test and solubility assay

- Measure OD<sub>600</sub> of each preculture
- Inoculate an 5 ml expression culture (LB+Kan [50 µg/ml]) with an OD<sub>600</sub> of 0.3
- Grow the expression culture for 1 hour at 37 °C shaking at 220 rpm.
- Measure OD<sub>600</sub> again (Optimal for induction of protein expression is an OD<sub>600</sub> between 0.4- 1.0)
- Take a 1 ml sample in an uninduced stage and harvest the cell by centrifugation at 8000 x g for 5 minutes
- Induce protein expression by adding IPTG (Final Concentration: 1mM).
- Incubate the culture for 3 hours at 37 °C shaking at 220 rpm.
- Take 1 ml sample 3 hours after induction and harvest the cells by centrifugation at 8000 x g for 5 minutes
- Harvest the rest 3 ml by centrifugation at 8000 x g for 5 minutes → used for solubility assay
- Keep the cell pellets at 4 °C
- Treat the negative control the same way!

### **25. August 2015**

#### 1) SDS-PAGE and Coomassie staining to verify successful pilotexpression of MEDH2


- Prepare Samples for SDS-Page:
  - Add 100 µl 1x SDS-Loading Buffer
  - Boil the sample at 95 °C for 20 minutes
- Run SDS-Page
  - 12 % Separating Gel
  - load 10 µl
- Stain for 1 hour in Coomassie Staining Solution
- Destain overnight



**Figure 1: Pilotexpression of MEDH2 in *E. coli* BL21.** Coomassie staining of heterologously expressed MEDH2 in pET-28. 1 ml samples were taken in an uninduced stage (n.i.) and 3 hours after induction (a.i.) with 1 mM IPTG. Theoretical weights: MEDH2- 40.15 kDa. Orange arrows determine the appropriate protein. As ladder PageRuler™ Plus Prestained Protein Ladder, 10 to 250 kDa (ThermoScientific) was used. As negative control (NC) selfligated pET-30 was used for expression in *E. coli* BL21. Proteins were separated by SDS-PAGE.

## 2) Solubility Assay, SDS-PAGE and Coomassie staining to verify soluble expression of MEDH2

- Solubility Assay
  - Resuspend the cell pellet in 600  $\mu$ l PBS buffer
  - Add glass beads
  - Use machine to disrupt the cells
  - Centrifuge at 16,000 x g for 5 minutes
  - Transfer 600  $\mu$ l supernatant to a new tube → soluble fraction
  - Resuspend the pellet in 300  $\mu$ l PBS buffer → Insoluble fraction
- Prepare Samples for SDS-Page:
  - Add 200  $\mu$ l 4x SDS-Loading Buffer
  - Boil the sample at 95 °C for 20 minutes
- Run SDS-Page
  - 12 % Separating Gel
  - load 10  $\mu$ l
- Stain for 1 hour in Coomassie Staining Solution
- Destain overnight



**Figure 2: Solubility assay of MEDH2.** Coomassie staining of solubility assay of MEDH2. Protein expression was performed for 3 hours. Cell pellet was lysed and separated into soluble (sol) and insoluble (insol) fraction. As control 1 ml sample was taken in an uninduced stage. Theoretical molecular weight: MEDH2-40.15 kDa. Orange arrows determine the appropriate protein. As ladder PageRuler™ Plus Prestained Protein Ladder, 10 to 250 kDa (ThermoScientific) was used. As negative control (NC) selfligated pET-28 was used for expression in *E. coli* BL21.