

I .Oligo

1. PCR System

Name	Total volume: 25 μ L	Total volume: 50 μ L
10*rTaq Buffer	2.5 μ L	5 μ L
2.5MdNTP	2 μ L	4 μ L
Oligo mix((1 μ M, Each primer 1/36 μ M)	9 μ L	18 μ L
rTaq	0.13 μ L	0.25 μ L
ddH2O	Add to total volume.	Add to total volume.

1. Procedure

94°C, 5min
94°C, 30s
55°C, 30s
72°C, 1min/kb
72°C, 7min
4-12°C, ∞

} $\times 30$

3. Gel Analysis

II .Overlap PCR Protocol 1

1. Amplify both two DNA fragments by conventional PCR.
2. Mix 50-100ng fragment I with equimolar fragment II in an Eppendorf tube.
3. Common PCR using the mixture in step 2 as templates.
4. Gel analysis.

II . Overlap PCR Protocol 2

1. Amplify both two DNA fragments by conventional PCR.
2. Mix 50-100ng fragment I with equimolar fragment II in an Eppendorf tube.
3. Conventional PCR without primers using the mixture in step 2 as templates.
Run for 8-12 cycles.
4. Conventional PCR using the product in step 3.Run for 25-35 cycles.
5. Gel analysis.