

Q5 site-directed mutagenesis (NEB)

Materials:

- Q5 hot start high-fidelity 2x master mix (12.5 μ l)
- 10 μ M forward primer (1.25 μ l)
- 10 μ M reverse primer (1.25 μ l)
- 1-25 ng μ l⁻¹ template DNA (1 μ l)
- 2x KLD reaction buffer (5 μ l)
- 10x KLD enzyme mix (1 μ l)
- Nuclease-free water

Protocol:

1. Assemble the following reaction in a PCR tube:

Q5 hot start high-fidelity 2x master mix	12.5 μ l
10 μ M forward primer	1.25 μ l
10 μ M reverse primer	1.25 μ l
Template DNA	1 μ l
Nuclease-free water	n μ l
Total volume	25 μ l

2. Mix reagents completely.
3. Perform the following cycling conditions:

Initial denaturation	98C	30 seconds
25 cycles	98C	10 seconds
	50-72C (Ta for primer set)	10-30 seconds
	72C	20-30 seconds per kb
Final extension	72C	2 mins
Hold	4-10C	Infinity

4. Use the PCR product to assemble the following reaction to allow kinase, ligase, and DpnI treatment:

PCR product	1 μ l
2x KLD reaction buffer	5 μ l
10x KLD enzyme mix	1 μ l
Nuclease-free water	3 μ l

5. Mix well with a pipette and incubate at room temperature for 5 mins. Use 5 μ l of the product to transform.