


Assembly:

PI_I13507



1st Day:

EXSP Digestion (see **Enzymatic Digestion Protocol**)

| | Part | Size | ng/μl |
|---|------------|---------|-------|
| 1 | PI_RBS | 212 bp | 104.8 |
| 2 | GFP_I13507 | 2585 bp | 273.5 |

| | Volume to 1,0 μg (μl) | Buffer 10x (μl) | BSA | Enzyme 1 | Volume (μl) | Enzyme 2 | Volume (μl) | H ₂ O to 20μl (μl) |
|---|-----------------------|-----------------|-----|----------|-------------|----------|-------------|-------------------------------|
| 1 | 9.5 | 2 (M) | - | SpeI | 1 | PstI | 1 | 10.5 |
| 2 | 4 | 2 (M) | - | XbaI | 1 | PstI | 1 | 12 |

| Final Plasmid | Resistance |
|---------------|------------|
| pSB1A2 | ampicillin |

Gel purification

- See PureLink® Quick Gel Extraction Kit Invitrogen™ manual
- Quantify digestion products

| Parts | ng/μl |
|------------|-------|
| PI_RBS | 9.5 |
| GFP_I13507 | 12.3 |

Obs: 260/280 is a quality parameter that tells you if your sample is contaminated with proteins. The greater it is compared to 1 the less contaminants you have.

Ligation (see **Ligation Protocol**)

| | | |
|-----------------------------|------------|--------------|
| Part containing the plasmid | Pl_RBS | 5 μ l |
| Insert | GFP_I13507 | 10.5 μ l |
| 10x T4 DNA Buffer | | 4 μ l |
| T4 DNA ligase 1u | | 0.5 μ l |
| H2O to 20 μ l | | - |

Obs: To determinate the amount of DNA necessary we used the following equation

$$\text{Insert ng} = \text{plasmid ng} \times \text{insert bp} / \text{plasmid bp} \times \text{insert: plasmid ratio}$$

- Incubate overnight at 37°C.
- Prepare and sterilize in the autoclave tubes with 6 ml of liquid LB medium.
- Prepare glycerol 40%

2nd Day:

Transformation (see **Transformation Protocol in Escherichia coli DH5- α**)

- Organism: E. coli DH5- α
- Selection: Ampicillin

4th Day:

Confirmation with NotI