


Assembly:

Pa_I13507



1st Day:

EXSP Digestion (see **Enzymatic Digestion Protocol**)

	Part	Size	ng/ μ l
1	Pa_RBS	161 bp	94.8
2	GFP_I13507	2585 bp	273.5

	Volume to 1,0 μ g (μ l)	Buffer 10x (μ l)	BSA	Enzyme 1	Volume (μ l)	Enzyme 2	Volume (μ l)	H ₂ O to 20 μ l (μ l)
1	10	2 (M)	-	SpeI	1	PstI	1	6
2	4	2 (M)	-	XbaI	1	PstI	1	12

Final Plasmid	Resistance
pSB1A2	ampicillin

Gel purification

- See PureLink® Quick Plasmide Miniprep Invitrogen™ manual
- Quantify digestion products

Parts	ng/ μ l
Pa_RBS	14
GFP_I13507	12.3

Obs: 260/280 in a quality parameter that tells you if your sample is contaminated with proteins. The greater it is compared to 1 the less contaminants you have.

Ligation (see **Ligation Protocol**)

Part containing the plasmid	Pa_RBS	4 μ l
Insert	GFP_I13507	10.5 μ l
10x T4 DNA Buffer		4 μ l
T4 DNA ligase 1u		0.5 μ l
H2O to 20 μ l		1

Obs: To determinate the amount of DNA necessary we used the following equation

$$\text{Insert ng} = \text{plasmid ng} \times \frac{\text{insert bp}}{\text{plasmid bp}} \times \text{insert: plasmid ratio}$$

- Incubate overnight at 37°C.
- Prepare and sterilize in the autoclave tubes with 6 ml of liquid LB medium.
- Prepare glycerol 40%

2nd Day:

Transformation (see **Transformation Protocol in Escherichia coli DH5- α**)

- Organism: E. coli DH5- α
- Selection: Ampicillin

4th Day:

Confirmation with NotI