


## Assembly:

PI\_I13507



## 1<sup>st</sup> Day:

### EXSP Digestion (see **Enzymatic Digestion Protocol**)

	Part	Size	ng/ $\mu$ l
1	PI_RBS	212 bp	104.8
2	GFP_I13507	2585 bp	273.5

	Volume to 1,0 $\mu$ g ( $\mu$ l)	Buffer 10x ( $\mu$ l)	BSA	Enzyme 1	Volume ( $\mu$ l)	Enzyme 2	Volume ( $\mu$ l)	H <sub>2</sub> O to 20 $\mu$ l ( $\mu$ l)
1	9.5	2 (M)	-	SpeI	1	PstI	1	10.5
2	4	2 (M)	-	XbaI	1	PstI	1	12

Final Plasmid	Resistance
pSB1A2	ampicillin

### Gel purification

- See PureLink® Quick Plasmide Miniprep Invitrogen™ manual
- Quantify digestion products

Parts	ng/ $\mu$ l
PI_RBS	9.5
GFP_I13507	12.3

**Obs:** 260/280 is a quality parameter that tells you if your sample is contaminated with proteins. The greater it is compared to 1 the less contaminants you have.

## Ligation (see **Ligation Protocol**)

Part containing the plasmid	Pl_RBS	5 $\mu$ l
Insert	GFP_I13507	10.5 $\mu$ l
10x T4 DNA Buffer		4 $\mu$ l
T4 DNA ligase 1u		0.5 $\mu$ l
H2O to 20 $\mu$ l		-

**Obs:** To determinate the amount of DNA necessary we used the following equation

$$\text{Insert ng} = \text{plasmid ng} \times \text{insert bp} / \text{plasmid bp} \times \text{insert: plasmid ratio}$$

- Incubate overnight at 37°C.
- Prepare and sterilize in the autoclave tubes with 6 ml of liquid LB medium.
- Prepare glycerol 40%

### **2<sup>nd</sup> Day:**

---

Transformation (see **Transformation Protocol in Escherichia coli DH5- $\alpha$** )

- Organism: E. coli DH5- $\alpha$
- Selection: Ampicillin

### **4<sup>th</sup> Day:**

---

Confirmation with NotI