


Assembly:

Pt_E0840



1st Day:

EXSP Digestion (see **Enzymatic Digestion Protocol**)

| | Part | Size | ng/ μ l |
|---|-------|--------|-------------|
| 1 | R0040 | 54 bp | 164.3 |
| 2 | E0840 | 878 bp | 171.0 |

| | Volume to 1,0 μ g (μ l) | Buffer 10x (μ l) | BSA (μ l) | Enzyme 1 | Volume (μ l) | Enzyme 2 | Volume (μ l) | H ₂ O to 50 μ l (μ l) |
|---|----------------------------------|-----------------------|----------------|----------|-------------------|----------|-------------------|---|
| 1 | 12.2 | 2 (M) | - | S | 1 | P | 1 | 3.80 |
| 2 | 11.7 | 2 (M) | 2 | X | 1 | P | 1 | 2.30 |

| Final Plasmid | Resistance |
|---------------|-----------------|
| pSB1C3 | chloramphenicol |

Gel purification

- See PureLink® Quick Plasmide Miniprep Invitrogen™ manual
- Quantify digestion products

| Parts | ng/ μ l |
|-------|-------------|
| R0040 | 5.3 |
| E0840 | 7.3 |

Obs: 260/280 in a quality parameter that tells you if your sample is contaminated with proteins. The greater it is compared to 1 the less contaminants you have.

Ligation (see **Ligation Protocol**)

| | | |
|-----------------------------|-------|-----|
| Part containing the plasmid | R0040 | 9.3 |
| Insert | E0840 | 8.5 |
| 10x T4 DNA Buffer | | 2 |
| T4 DNA ligase 1u | | 0.4 |
| H2O to 20µl | | - |

Obs: To determinate the amount of DNA necessary we used the following equation

$$\text{Insert ng} = \text{plasmid ng} \times \frac{\text{insert bp}}{\text{plasmid bp}} \times \text{insert: plasmid ratio}$$

- Incubate overnight at 37°C.
- Prepare and sterilize in the autoclave tubes with 6 ml of liquid LB medium.
- Prepare glycerol 40%

2nd Day:

Transformation (see **Transformation Protocol in Escherichia coli DH5-α**)

- Organism: E. coli DH5-α
- Selection: Ampicillin

4th Day:

Confirmation with NotI.