



UNSW
AUSTRALIA

BABS UNSW iGEM Lab Protocol



Procedure	Name		Mammalian cell invasion assay			
	Description		Testing of the IPTG induced invasion+LLO+GFP system across all mammalian and bacterial cells.			
Document	Name	Kristen Perry	Date	20/08/15	Version	1
Requirements	Time	Induction and incubation – 2 days Spec. prep. – 30 mins + 24 hr incubation Microscopy – 1 day				
	PPE	Gloves, lab coats (experiments) Gloves, gowns (microscopy)				
	Equipment	210 x small culture dishes 210 x 0.17 mm round coverslips 105 x microscope slides Pipettes, tips Centrifuge OD reader				
	Materials	Mammalian cells Bacterial cells transformed with IPTG+INV+LLO+GFP L-15 media L-15 media + 100 µg/ml gentamicin IPTG PBS Fixing and mounting materials – 0.3 M glycine, 3-4% paraformaldehyde, Prolon gold				
Step 1	Prepare 210 small cell culture dishes with mammalian cells in L-15 media, with coverslips inserted in the base I know this is a lot, and that's just for one mammalian cell type – this would be ideal in my opinion, but if we need to cut down, we could remove the 10:1 and/or 8:1 ratio, or the 18 hr time point					
Step 2	Measure concentration/confluence of mammalian cells and calculate cell population					
Step 3	Pellet transformed bacteria, gently wash in PBS and resuspend					

Step 4	Dilute to desired OD (measured per 100 ul) such that ratios range from <1:1 to 10:1																																																																																																																																																																																																																		
Step 5	Gently wash mammalian cells in L-15 media (and resuspend (EDTA/trypsin)?) in L-15 Unsure as to whether we should resuspend monolayer to provide more cell surface area at this step																																																																																																																																																																																																																		
Step 6	Add 100 uL of bacterial solution to each well (100 uL PBS for control) and add IPTG to induce gene expression																																																																																																																																																																																																																		
Step 7	<p>Incubate at 37 degrees in atmospheric CO2 as follows:</p> <p>Synthetic chloroplast paper incubated their cyanobacteria cultures at 30 degrees, unsure whether this is a good idea? Optimal temp. seems to be 30-33 for Synechocystis</p> <p>e.g. HeLa cells</p> <table border="1" data-bbox="464 909 1434 1021"> <tr><td>E. coli</td></tr> <tr><td>Lactococcus</td></tr> <tr><td>Synechocystis</td></tr> </table> <table border="1" data-bbox="464 1061 1434 2018"> <thead> <tr> <th colspan="2"></th> <th colspan="7">RATIO OF BACTERIA TO HOST</th> </tr> <tr> <th rowspan="2">INCUBATION TIME</th> <th rowspan="2">3 hour incubation</th> <th>0</th> <th><1:1</th> <th>2:1</th> <th>4:1</th> <th>6:1</th> <th>8:1</th> <th>10:1</th> </tr> </thead> <tbody> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr> <th rowspan="2">6 hour incubation</th> <th>0</th> <th><1</th> <th>2</th> <th>4</th> <th>6</th> <th>8</th> <th>10</th> </tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr> <th rowspan="2">9 hour incubation</th> <th>0</th> <th><1</th> <th>2</th> <th>4</th> <th>6</th> <th>8</th> <th>10</th> </tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr> <th rowspan="2">18 hour incubation</th> <th>0</th> <th><1</th> <th>2</th> <th>4</th> <th>6</th> <th>8</th> <th>10</th> </tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr><td>0</td><td><1</td><td>2</td><td>4</td><td>6</td><td>8</td><td>10</td></tr> <tr> <th></th> <th></th> <th>0</th> <th><1</th> <th>2</th> <th>4</th> <th>6</th> <th>8</th> <th>10</th> </tr> <tr> <td></td> <td></td> <td>0</td> <td><1</td> <td>2</td> <td>4</td> <td>6</td> <td>8</td> <td>10</td> </tr> </tbody> </table>	E. coli	Lactococcus	Synechocystis			RATIO OF BACTERIA TO HOST							INCUBATION TIME	3 hour incubation	0	<1:1	2:1	4:1	6:1	8:1	10:1	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	6 hour incubation	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	9 hour incubation	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	18 hour incubation	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10	0	<1	2	4	6	8	10			0	<1	2	4	6	8	10			0	<1	2	4	6	8	10
E. coli																																																																																																																																																																																																																			
Lactococcus																																																																																																																																																																																																																			
Synechocystis																																																																																																																																																																																																																			
		RATIO OF BACTERIA TO HOST																																																																																																																																																																																																																	
INCUBATION TIME	3 hour incubation	0	<1:1	2:1	4:1	6:1	8:1	10:1																																																																																																																																																																																																											
		0	<1	2	4	6	8	10																																																																																																																																																																																																											
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
6 hour incubation	0	<1	2	4	6	8	10																																																																																																																																																																																																												
	0	<1	2	4	6	8	10																																																																																																																																																																																																												
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
9 hour incubation	0	<1	2	4	6	8	10																																																																																																																																																																																																												
	0	<1	2	4	6	8	10																																																																																																																																																																																																												
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
18 hour incubation	0	<1	2	4	6	8	10																																																																																																																																																																																																												
	0	<1	2	4	6	8	10																																																																																																																																																																																																												
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
0	<1	2	4	6	8	10																																																																																																																																																																																																													
		0	<1	2	4	6	8	10																																																																																																																																																																																																											
		0	<1	2	4	6	8	10																																																																																																																																																																																																											

		24 hour incubation	0	<1	2	4	6	8	10
			0	<1	2	4	6	8	10
			0	<1	2	4	6	8	10
			0	<1	2	4	6	8	10
			0	<1	2	4	6	8	10
Step 8	Following incubation, wash cells three times with PBS and replace media with L-15 containing 100 ug/mL gentamycin (does not permeate mammalian cell membrane)								
Step 9	Rinse cells in PBS and remove supernatant								
Step 10	Fix immediately in 3-4% paraformaldehyde for 10-20 minutes, with 0.3M glycine added. Rinse briefly with PBS.								
Step 11	Mount with Prolon gold and leave for 24 hours.								
Step 13	Image with brightfield Zeiss microscope. Before imaging, dab coverslip clean of PBS with Kimiwipe and water to avoid crystal artefact								
Notes	Need to confirm whether IPTG in cultures will autofluoresce within GFP emission or otherwise affect imaging								
Version History	First version – up for correction and annotations								