

# P18: FORMATION OF LIPOSOMES

## Aim: design of photopolymerizable phospholipids

### Reagents

- ✦ Photopolymerizable phospholipids
- ✦ 140mM NaCl
- ✦ 10mM HEPES Buffer

### Materials

- ✦ Pipettes
- ✦ Pipettes tips
- ✦ Centrifuge
- ✦ Vortex
- ✦ Centrifuge
- ✦ Vacuum desiccator

### Protocol

#### **This is a very light protocol we made from reference article**

- ✦ Mix lipids at desired molar ratios in a glass tube.
- ✦ Remove the solvent under nitrogen to form a lipid film.
- ✦ Remove residual chloroform by placing the films overnight in a vacuum desiccator.
- ✦ Formation of multilamellar vesicles by reconstituting the lipid with HEPES buffer (10mM HEPES, 140mM NaCl, pH 7.5) and by vigorous vortexing.
- ✦ To encapsulate toxins, reconstitute lipid films with the desired concentrations of toxins in HEPES buffer.
- ✦ Sonication at 4 degrees for 5-10 minutes (1 minute pulse and 1 minute rest)
- ✦ Centrifuge the samples at 2000g for 5-10 minutes to remove any titanium particles and larger aggregates.
- ✦

Je mets quoi d'autre...?