# LB Medium

## Liquid medium (Total volume: 1000mL)

1L -Water -Peptone 10g -NaC1 10g -Yeast extract powder

# Solid medium (Total volume: 1000mL)

1L -Water -Peptone 10g -NaC1 10g -Yeast extract powder 5g -Agar powder 15g

### LB Minimal Salts Medium

# Liquid medium (Total volume: 1000mL)

-Water 1L -Peptone 10g

-NaCl 5g

-Yeast extract powder 5g

-pH should be 7.0

#### medium(Total volume:1000mL) Solid

-Water 1L -Peptone 10g -NaC1 5g -Yeast extract powder 5g -Agar powder 15g

-pH should be 7.0